

A close-up photograph of an elderly woman with short, light-colored hair lying in a hospital bed. She is looking towards the left of the frame. A dark-colored blood pressure cuff is wrapped around her left upper arm. The cuff has white text on it. The background is slightly blurred, showing a person in a white coat, likely a healthcare professional, standing by the bed. The overall lighting is soft and clinical.

Durable One-Piece Blood Pressure Cuffs

WelchAlllyn®

Advancing Frontline Care™

WelchAlllyn®
ADULT II
DURABLE
Blood Pressure Cuff



Supreme durability. Simple design.

Welch Allyn continues to set the standard in quality and durability with our extensive line of Durable One-Piece Blood Pressure Cuffs. Long-lasting, easy to use, and affordable, they're the cuff of choice for hospitals and clinics around the world. That's because we offer virtually every type of cuff an institution may need, along with the ability to simplify inventory like no other cuff provider.

Our Durable One-Piece Cuffs meet the latest clinical guidelines for proper fit from AAMI and AHA, and because we manufacture the cuffs ourselves using a fully automated process, they are always of the highest quality.



A level of quality only Welch Allyn can offer.

For nearly a century, Welch Allyn has made blood pressure instruments with unmatched quality. So when we began manufacturing Durable One-Piece Blood Pressure Cuffs, we adhered to the same standards healthcare providers have come to expect.

- Made with soft but durable 70-denier nylon
- Tested to last thousands of cycles; 3-year warranty
- Validated for use on multiple monitors
- Patented flat port tubing eliminates leaks
- No sewn seams to tear or fray

A safe, smart choice for all healthcare facilities.

With no separate bladder to remove, clean, or replace, Welch Allyn Durable One-Piece Blood Pressure Cuffs deliver safe and effective performance throughout their long life cycle.

- Easy to clean with soap and water; disinfect with glutaraldehyde, or alcohol or bleach solutions; and gas sterilize
- Antimicrobial treated to help combat mildew and bacteria growth throughout the life of the cuff
- 100% latex-free cuffs minimize risk of allergic reaction

A streamlined inventory for simplicity's sake.

Because Welch Allyn offers such an extensive line of blood pressure cuffs, you can make all your purchases from a single source. And Welch Allyn sales representatives can help you decide how best to reduce your cuff inventory. That could mean providing cuffs that fit any brand of monitors your institution has, or standardizing all monitors throughout the facility with one tubing and connector configuration.

Either way, receiving all cuffs from a single supplier means simplified inventory, fewer purchase orders, and lower costs. And it means reducing the time and frustration practitioners face locating the appropriate cuff.


- Every cuff size needed – from newborn to thigh
- One- and two-tube varieties
- Connectors attached to fit virtually any electronic NIBP monitor or sphygmomanometer
- Ability to standardize with one tubing and connector configuration

Welch Allyn Durable One-Piece Blood Pressure Cuffs can be customized to work with almost any type of electronic vital signs monitor.




Ability to fit any electronic NIBP monitor or sphygmomanometer.

Welch Allyn cuffs are shipped with connectors attached that allow them to be used with virtually any NIBP monitor or sphygmomanometer. Durable One-Piece Cuffs are available with any of the following connectors in either one- or two-tube versions.




Male Screw w/ Female Luer Slip

Fits Welch Allyn Atlas, Vital Signs Monitors, and Spot (1 tube) and Critikon Dinamap Monitors (2 tubes)



Female Luer Lock

Fits Sphygmomanometers and BCI, Colin, and Datascope Monitors (1 tube)



Male Bayonet

Fits HP (Philips), Spacelabs, Datascope, Siemens, and Criticare Monitors (1 tube)



Female Locking

Fits Welch Allyn Protocol (1 tube) and GE Marquette Monitors (2 tubes)

one-tube

	Screw (each)	Screw (5 pack)	Luer Lock (each)	Luer Lock (5 pack)	Bayonet (5 pack)	Locking (5 pack)
• NEWBORN (7.7-10.5 CM)	5082-201-3	5082-81-3	5082-201-1	5082-81-1	5082-81-3HP	5082-81-3MQ
• INFANT (9.8-13.3 CM)	5082-202-3	5082-82-3	5082-202-1	5082-82-1	5082-82-3HP	5082-82-3MQ
• SMALL CHILD (12.4-16.8 CM)	5082-203-3	5082-83-3	5082-203-1	5082-83-1	5082-83-3HP	5082-83-3MQ
• CHILD (15.8-21.3 CM)	5082-204-3	5082-84-3	5082-204-1	5082-84-1	5082-84-3HP	5082-84-3MQ
• CHILD PRINT (15.8-21.3 CM)	5082-204P-3	5082-84P-3	5082-204P-1	5082-84P-1	5082-84P-3HP	5082-84P-3MQ
• SMALL ADULT (20.0-27.0 CM)	5082-205-3	5082-85-3	5082-205-1	5082-85-1	5082-85-3HP	5082-85-3MQ
• ADULT (25.3-34.4 CM)	5082-206-3	5082-86-3	5082-206-1	5082-86-1	5082-86-3HP	5082-86-3MQ
• ADULT LONG (25.3-34.3 CM)	5082-206L-3	5082-86L-3	5082-206L-1	5082-86L-1	5082-86L-3HP	5082-86L-3MQ
• LARGE ADULT (32.1-43.4 CM)	5082-207-3	5082-87-3	5082-207-1	5082-87-1	5082-87-3HP	5082-87-3MQ
• LG. ADULT LONG (32.1-43.4 CM)	5082-207L-3	5082-87L-3	5082-207L-1	5082-87L-1	5082-87L-3HP	5082-87L-3MQ
• THIGH (40.7-55.0 CM)	5082-208-3	5082-88-3	5082-208-1	5082-88-1	5082-88-3HP	5082-88-3MQ
MULTI-PACK (1 of each size, except Print and Longs)	5082-241-3	N/A	5082-241-1	N/A	5082-241-3HP	5082-241-3MQ

two-tube

	Screw (each)	Screw (5 pack)	Luer Slip* (each)	Luer Slip* (5 pack)	Bayonet (5 pack)	Locking (5 pack)	Inflation System (each)
• INFANT (9.8-13.3 CM)	5082-202-4	5082-82-4	5082-202-2	5082-82-2	5082-82-4HP	5082-82-4MQ	5082-222-1
• SMALL CHILD (12.4-16.8 CM)	5082-203-4	5082-83-4	5082-203-2	5082-83-2	5082-83-4HP	5082-83-4MQ	5082-223-1
• CHILD (15.8-21.3 CM)	5082-204-4	5082-84-4	5082-204-2	5082-84-2	5082-84-4HP	5082-84-4MQ	5082-224-1
• CHILD PRINT (15.8-21.3 CM)	5082-204P-4	5082-84P-4	5082-204P-2	5082-84P-2	5082-84P-4HP	5082-84P-4MQ	5082-224P-1
• SMALL ADULT (20.0-27.0 CM)	5082-205-4	5082-85-4	5082-205-2	5082-85-2	5082-85-4HP	5082-85-4MQ	5082-225-1
• ADULT (25.3-34.4 CM)	5082-206-4	5082-86-4	5082-206-2	5082-86-2	5082-86-4HP	5082-86-4MQ	5082-226-1
• ADULT LONG (25.3-34.3 CM)	5082-206L-4	5082-86L-4	5082-206L-2	5082-86L-2	5082-86L-4HP	5082-86L-4MQ	5082-226L-1
• LARGE ADULT (32.1-43.4 CM)	5082-207-4	5082-87-4	5082-207-2	5082-87-2	5082-87-4HP	5082-87-4MQ	5082-227-1
• LG. ADULT LONG (32.1-43.4 CM)	5082-207L-4	5082-87L-4	5082-207L-2	5082-87L-2	5082-87L-4HP	5082-87L-4MQ	5082-227L-1
• THIGH (40.7-55.0 CM)	5082-208-4	5082-88-4	5082-208-2	5082-88-2	5082-88-4HP	5082-88-4MQ	5082-228-1
MULTI-PACK (1 of each size, except Print and Longs)	5082-241-4	N/A	5082-241-2	N/A	5082-241-4HP	5082-241-4MQ	N/A

*For European-style cuffs with two long tubes, replace "-2" with "-5" in the part number (e.g., 5082-202-5), excluding 5082-241-2.

