


Versatile Illumination Systems for Optimal Visualisation



Sigmoidoscope/Anoscope Lighting Systems

Powerful illumination systems.

- > Halogen light provides consistent, long-lasting illumination (73500, 73550 only)
- > Fibre optics provide a cool distal ring of light with no reflections or obstructions (36019, 38700 only)
- > Sigmoidoscopes/anoscopes can be used with a choice of illumination systems - a 6 V system with optional on/off switch and rheostat control, or with Welch Allyn's 3.5 V/2.8 V or 2.5 V hand-held power sources (with adaptors 73500 and 73550).
- > 73500 adapter fits most Welch Allyn 3.5 V power handles
- > 73550 adapter fits most Welch Allyn 2.5 V power handles (except PocketScope™)


Sigmoidoscope/Anoscope Lighting Systems and Accessories

30200	Insufflation Bulb, Complete	73210	Light Handle with 0.9m cord
30130	30 cm Suction Tube (Standard)	73308	Optional On/Off Switch
30140	40 cm Suction Tube	73324	Transformer with 1.5m cord
36019	Fibre Optic Light Head for Disposable Sigmoidoscopes and Anoscopes	73400	Optional Rheostat Control
36146	Illumination System for Disposable Sigmoidoscopes (Includes: 73324, 73210, 36019 and 30200)	73500	3.5 V Halogen Handle Adapter
38300	2.5 V Vacuum Light Carrier	73550	2.5 V Halogen Handle Adapter
38700	Fibre Optic Light Carrier for use with Standard Anoscopes and Light Adapter	00200-U	2.5 V Vacuum Lamp (small) (for 38300)
71020-A	3.5 V Rechargeable Handle	03100-U	3.5 V Halogen Lamp (for 73500)
		03400-U	2.5 V Halogen Lamp (for 73550)
		07800-U	6 V Halogen Lamp (for 73210)

For more information, contact us at 0207 365 6780 today!

Welch Allyn UK Ltd,
Cublington Road,
Aston Abbots,
Buckinghamshire HP22 4ND
Tel: 0207 365 6780 Fax: 0207 365 9694 www.welchallyn.com

WelchAllyn
Advancing Frontline Care™



RIGID ENDOSCOPY

WelchAllyn
Advancing Frontline Care™

Quality and Value in Rigid Endoscopy

Welch Allyn Disposable Sigmoidoscopes

Welch Allyn offers a wide range of rigid endoscopy products for anoscopy and sigmoidoscopy, complemented by versatile fibre optic illumination systems for optimum visualisation.

Welch Allyn disposable sigmoidoscopes are available in a variety of configurations to accommodate all your clinical needs and come complete with obturator making them ready for immediate use. Illumination is provided by a 6 V illumination system.

- 53130** Disposable Sigmoidoscopes (4 boxes of 25)
- 53130-L** Lubricated Disposable Sigmoidoscopes (4 boxes of 25)
- 53130-F** Disposable Sigmoidoscopes (4 boxes of 25) each with filter
- 53130-LF** Lubricated Disposable Sigmoidoscopes (4 boxes of 25) each with filter
- 36146** Complete Illumination System for Disposable Sigmoidoscope comprising fibre optic light head, light handle, transformer and insufflation bulb

Filters may also be purchased separately (see next page for details).



Welch Allyn also offers a wide range of reusable standard and fibre optic sigmoidoscopes—please contact us for further information on sizes and product codes.



Insufflation Systems and Filters

- > To aid clinical diagnosis and help prevent cross contamination, the Welch Allyn sigmoidoscopes should be used in conjunction with the insufflation system and filter.
- > Insufflation helps distend the examination area to aid viewing
- > Filter prevents transmission of matter between patients
- > Filters to be used in conjunction with Luer Lock adaptor (available separately).

- 30200** Insufflation Bulb, complete
- 30210** Hydrophobic filters (pack of ten)
- 30215** Luer Lock for use with hydrophobic filters (pack of ten)