

See It. Share It. Save It.



Welch Allyn Digital MacroView™ Otoscope

The Welch Allyn Digital MacroView Otoscope is a simple, affordable tool that lets you see, explain, and save digital images of the ear.

Offering you clear, live and still images of the ear canal and tympanic membrane, the Digital MacroView™ Otoscope is a valuable tool in support of the treatment and management of your patients' ears. Easy to learn and use, the Digital MacroView is a plug and play system requiring no additional purchase of software.

Here's how the Welch Allyn Digital MacroView Otoscope can benefit you:

- > Improves communication which can lead to better patient compliance
- > Allows you to store digital images in patient records for enhanced documentation of exam findings
- > Provides you the ability to share images with other providers as visual descriptive reference
- > Gives you an excellent training aid for educational purposes
- > Offers you the ability to view small details of tubes and hearing aids



*handle sold separately.
(Laptop computer not included).

WelchAllyn®

Advancing Frontline Care™

The Welch Allyn Digital MacroView™ Otoscope

Features you expect from an otoscope and more



USING THE DIGITAL MACROVIEW OTOSCOPE IS EASY AS 1 – 2 – 3

1 Plug In.

Plug the otoscope into your computer using the USB cable.

2 Power On.

Put on ear specula and turn on light.

3 Play.

View live picture or capture image.

Ordering Information

23920 Digital MacroView Otoscope with 3 m USB cable and set of four reusable ear specula

Accessories

23804	Insufflator Bulb	52434-U	4.25 mm Single use Ear Specula,
24400-U	Additional set of four Reusable Specula	06500-U	Lamp
24302-U	2.5 mm Universal Poly Reusable Speculum	23900-10	Replacement USB Cable, 3 m
24303-U	3.0 mm Universal Poly Reusable Speculum	71062-A	NiCad Rechargeable Handle and Transformer - Europe
23404-U	4.0 mm Universal Poly Reusable Speculum	71064-A	NiCad Rechargeable Handle and Transformer - UK
23405-U	5.0 mm Universal Poly Reusable Speculum	71062-C	NiCad/Convertible Rechargeable Handle and Transformer - Europe
24320	SofSeal Small, Disposable Sealing Tip	71064-C	NiCad/Convertible Rechargeable Handle and Transformer - UK
24330	SofSeal Medium, Disposable Sealing Tip	71902	Lithium Ion Handle, Charger and Transformer - Europe
52432-U	2.75 mm Single use Ear Specula	71904	Lithium Ion Handle, Charger and Transformer - UK

System Requirements: Windows® XP Service Package 2 or higher with an available USB 2.0 port.

To find out how the new Welch Allyn Digital MacroView Otoscope can help you see, explain and save digital images of the ear, go to www.welchallyn.com or call +353 46 90 67790.

Welch Allyn UK Ltd, Cublington Road, Aston Abbots, Buckinghamshire HP22 4ND, England
Tel: +44 (0) 207 365 6780 Fax: +44 (0) 207 365 9694

Welch Allyn South Africa (Pty) Ltd, 394 Surrey Avenue, Ferndale, Randburg, 2194, South Africa
P O Box 809, Fountainsbleau, South Africa 2032
Tel: +27 (0) 11 777-7555 Fax: +27 (0) 11 777-7556

Welch Allyn Customer Services, Navan Business Park, Dublin Road, Navan, Co. Meath, Republic of Ireland
Tel: +353 46 90 67790 Fax: +353 46 90 67754

WelchAllyn®

Advancing Frontline Care™