



# Welch Allyn Propaq® CS and Propaq® LT Vital Signs Monitors

Cost-effective monitoring solutions for outpatient surgery and procedures.

**WelchAllyn®**

Advancing Frontline Care™

# Easy-to-use. Reliable. Cost-effective.

And versatile enough to be used anywhere.

Welch Allyn Propaq CS and Propaq LT Patient Monitors

## Easy-to-use

Designed with the busy user in mind, the Propaq CS features a bright, full-color touchscreen display and easy-to-use menu functions, making it an ideal monitor during procedures requiring conscious sedation and continuous patient monitoring.

Ideal for continuous monitoring and spot-check during simple procedures, the Propaq LT features an intuitive, navigational keypad with the option to connect to a large color display for across-the-room viewing of patient data. Plus, the Propaq LT weighs in at just about two pounds, making it easy to move throughout your facility.

## Reliable enough for any situation

The Propaq product line was designed for rugged military use and can stand up to your toughest days. The Propaq LT has shock resistance that can withstand a six-foot drop, so you can be confident that it is durable enough to use during simple procedures or monitoring vitals in pre-op and recovery. And with legendary Propaq quality, you can rely on your Propaq LT and Propaq CS monitors for many years to come.

## Cost-effective solutions

When you purchase a Welch Allyn Propaq monitor, you're investing in a quality brand that has stood the test of time. For decades, both hospital and military personnel have counted on Propaq monitors for their proven reliability and durability, so you can be confident in your investment. And to further protect your investment, the Welch Allyn Partners in Care Services<sup>SM</sup> program will help ensure that your monitors are up and running whenever you need them.



## Welch Allyn Propaq® CS Patient Monitor

Recommended for continuous vital signs monitoring for procedural sedation and pre- and post-operative care.



Ideal for outpatient surgery and recovery:

Gastroenterology  
Endoscopy  
Pain Management  
Ophthalmology  
Oral/Maxillofacial  
Multispecialty clinics



Complete wall, pole, and roll stand mounting options and accessory storage baskets can be added to optimize the use of your Propaq CS Monitor.

Propaq CS features include:

3- or 5-lead ECG, heart rate
Masimo SET® or Nellcor® OxiMax® pulse oximetry
NIBP with Welch Allyn FlexiPort® technology
Pulse rate
Respiration
Temperature
Capnography (sidestream and/or mainstream)*
Invasive blood pressure (one or two channels)*
Full patient alarms
Apnea alarming
Trending on all parameters
Programmable defaults on displayed parameters and alarm limits
Printer*
RS 423 connectivity port

\*Optional parameters

## Welch Allyn Propaq® LT Patient Monitor



Ideal for spot check and continuous monitoring needs during simple procedures and pre- and post-operative care.



Print an 8.5" x 11" tabular trend report to show patient data, highlight alarm events, or create a PDF to attach to your electronic medical records from your PC.



Propaq LT is versatile enough to be used anywhere you need it most.


Propaq LT features include:

3- or 5-lead ECG, heart rate
Masimo SET® or Nellcor® OxiMax® pulse oximetry
NIBP with Welch Allyn FlexiPort® technology
Pulse rate
Respiration
Continuous monitoring and spot check modes
Storage of up to 20 full-disclosure data snapshots
Print 8.5" x 11" reports from your PC via USB
Ability to attach PDF to an electronic medical record from your PC
24-hour tabular and graphical trending
Compact design weighing about two pounds, around the size of a paperback book
Large color display option



## Welch Allyn legendary Propaq® quality

Welch Allyn Propaq® monitors have been used in hospitals and by the military for decades due to their proven quality, durability and ease-of-use in demanding and unpredictable situations. In fact, the Propaq LT has undergone rigorous testing and is certified for use during all phases of flight aboard U.S. Army aircraft.



## Partners in Care Services<sup>SM</sup>

SUPPORT | EDUCATION | PROFESSIONAL

### Get the Welch Allyn Partners in Care Services<sup>SM</sup> advantage.

Partners in Care goes beyond your standard product warranty to take service and support to another level in order to meet your goals for cost and life cycle planning.

#### Welch Allyn Partners in Care Services provides you with:

- Improved system uptime and performance
- Predictable support costs
- Regulatory compliance
- Faster response times and issue resolution
- Efficient, cost-effective methods to keep equipment at appropriate validation levels

### Comprehensive Partnership Program

When you have limited internal resources or simply want your staff to concentrate on other tasks without relying on third-party providers, consider our Comprehensive Partnership Program—it's like adding Welch Allyn service technicians to your staff.

#### We partner with you to handle all of your service needs including:

- Expedited service repair and delivery
- Free shipping
- Calibration discounts
- Software updates
- Loaner coverage\*
- Parts and labor at no charge
- Accessory protection

*\*Loaners are included under warranty with Propaq CS and Propaq LT patient monitors*

*Please discuss your needs with your Welch Allyn representative, and they can customize a service and support program that works best for your facility.*

Welch Allyn Protocol, Inc.  
8500 S.W. Creekside Place, Beaverton, OR 97008 USA  
(p) 503.530.7500 (f) 503.526.4200  
Technical Support: 800.289.2500  
[www.welchallyn.com](http://www.welchallyn.com)

**WelchAllyn**<sup>®</sup>

Advancing Frontline Care<sup>™</sup>

