

Approved PDAs for use with Welch Allyn[®] AcuityLink[®] Clinician Notifier

This list indicates the Welch Allyn approved PDA devices. If your device is not on this list, contact your Welch Allyn account representative.

Installation of Clinician Notifier software on a mobile device not specified by Welch Allyn can interfere with effective Clinician Notifier program performance, including alarm notification. Only use devices approved by Welch Allyn.

Approved PDAs

Motorola MC55/Symbol MC5590

- Operating system: Windows Mobile 6.1 classic
- Screen resolution: 320 x 240

Motorola MC50/MC5040 (upgrade only)

- Operating system: Windows mobile 2003 SE 4.21.1088
- Screen resolution: 320 x 240
- For information on upgrading Clinician Notifier on this PDA, see the instructions under the Technical tab at <http://www.welchallyn.com/apps/products/product.jsp?id=16-vo-96-1232546273407>, and click on the device name.

For PDAs listed here:

- **HHP 7900**
- **HHP 9500**
- **Honeywell Dolphin 7600**
- **Honeywell Dolphin 9900**
- **Motorola MC7090**
- **HP iPaq 111**
- **Pharos 600e**
- **Unitech PA600 MCA**

Use the following settings:

- Operating system: Pocket PC 2003 SE, Windows Mobile 5.x, Windows Mobile 6.0x or 6.1x (see device documentation)
- Screen resolution: 320 x 240
- For information on installing Clinician Notifier on these devices, see the instructions under the Technical tab at <http://www.welchallyn.com/apps/products/product.jsp?id=16-vo-96-1232546273407> and click on the device name.