

Welch Allyn SureTemp® Plus 690 Electronic Thermometer

FOR VETERINARIANS



Now taking temperatures is easier, safer and more accurate.

The Welch Allyn SureTemp Plus 690 is one of the fastest, most reliable and accurate thermometers ever developed for small and large animals alike, providing you with a 10- to 13-second rectal temperature measurement. Simple, yet sophisticated, it has an easy-to-read LCD display, leaving your hands free to work with your patients, and the interchangeable, removable probe well reduces the risk of cross-contamination.

WelchAllyn®

Advancing Frontline Care™

The SureTemp® Plus 690 Electronic Thermometer

With the speed and accuracy you require, and the versatility you need for use with all your patients, all the time.



- Interchangeable probe with push-button probe cover ejection

- C/F button for toggling between Celsius and Fahrenheit

- Prominent recall button for display of last temperature taken



- Convenient storage for 25 probe covers

- Removable probe well for easy cleaning

- Mode selection button for easy toggling between customization options in the biomed mode


- Ergonomic design with rubberized grips



- Optional wall mount has storage space for 50 additional probe covers



- Interchangeable rectal probe and well reduce risk of cross-contamination



- Large, easy-to-read LCD

- Battery Life indicator

Welch Allyn SureTemp Plus 690 Specifications and Ordering Information

Product Description	SURETEMP PLUS 690 THERMOMETER
Product Type	Electronic Thermometer
Calibration Accuracy	±0.2°F (±0.1°C) Monitor Mode
Temperature Predict Time	Approximately 10 to 13 seconds (Rectal)
Unit Dimensions	8.46" x 3.18" x 2.43" (21 mm x 81 mm x 62 mm)
Weight	12.6 ounces (357 grams)
Ambient Temperature Range	50°F – 104°F (10°C – 40°C)
Patient Temperature Range	80°F – 110°F (26.7°C – 43.3°C)
Power Source	(3) 1.5 VDC AA batteries
Battery Life	Approx. 6,000 temperatures
Warranty	2-year limited instrument 1-year limited probe 90-day limited probe well

01690-201V SureTemp Plus 690 Electronic Thermometer with Rectal Probe with 4 ft cord and Rectal Probe Well, 2-year limited warranty

01690-301V SureTemp Plus 690 Electronic Thermometer with Wall Mount, Rectal Probe with 9 ft cord and Rectal Probe Well, 2-year limited warranty

Accessories

- 21326-0000 Wall Holder for 690/692
- 21330-0000 Holder for Spare Probe and Probe Well for 690/692
- 06138-000 Calibration Key, Assembly for 690/692
- 02892-003 Probe and Well Kit, 4 ft cord on rectal probe
- 02892-103 Probe and Well Kit, 9 ft cord on rectal probe
- 02891-1000 Rectal Probe Well, red for 690/692
- 08273-000 Thermometer Stand with Holder for Spare Probe and Well
- 21333-0000 Bracket, M600 Stand Upgrade, with Holder for Spare Probe and Well
- 98022-000 Extended Warranty M690

Replacement Probe Covers

- 05031-101 Disposable Probe Covers (case of 1,000; packaged 25 per box)

For more information or to order, please contact your local Welch Allyn veterinary distributor, go to www.welchallyn.com/dvm or call 1.800.535.6663.

4341 State Street Road, Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361 www.welchallyn.com/dvm

WelchAllyn®

Advancing Frontline Care™