



Solid State Headlights From Welch Allyn

High-quality, lightweight solid state headlights made affordable for every veterinary practice



Solid State Procedure Headlight

- 50 lumens of bright, white, shadow-free light for true tissue color
- Cool beam, safe to touch
- At just six ounces, it's one of the most lightweight and comfortable headlights available
- Estimated life of up to 10,000 hours of continuous burn time
- Variable spot size for well-defined illumination
- Available with direct power source or the rechargeable "belt-clip" portable power source



Solid State Portable Headlight

- 20 lumens of cool, white light for true tissue color
- Weighing just four ounces, it's one of the most lightweight headlights available
- Estimated life of up to 10,000 hours of continuous burn time
- Variable spot size for well-defined illumination
- Rechargeable battery housed in soft, terry cloth headband allows the user freedom of movement

WelchAllyn[®]

Advancing Frontline Care[™]



Solid State Headlight Specifications



Solid State Procedure Headlight

Specifications

Light Source: 5-watt solid state lamp, color temperature 5500° Kelvin, 50 lumens, estimated life of up to 10,000 hours continuous burn

Working Distance: 16 inches

Spot Size Range: 1.5 to 6 inches at the working distance

Weight: 6 ounces

Battery Run Time: 1-hour continuous runtime on a fully charged battery

Battery Life: Approximately 2.5 years (depending on usage)

Warranty: 1-year limited warranty

Model Numbers

49020 Solid State Procedure Headlight with Direct Power Supply

Accessory Numbers

49099 Carrying Case

75260 Battery Pack with Charger

Replacement Parts

49000 Procedure Headlight only

72250 Replacement Battery Pack

74180 Charger

75270 Direct Power Supply (U.S.)

Solid State Portable Headlight

Specifications

Light Source: 1-watt solid state lamp, color temperature 5500° Kelvin, 20 lumens, estimated life of up to 10,000 hours continuous burn

Working Distance: 16 inches

Spot Size Range: 1.5 to 6 inches at the working distance

Weight: 4 ounces

Battery Run Time: Approximately 90 minutes of continuous run-time on a fully charged battery

Battery Life: Approximately 2 years (depending on usage)

Warranty: 1-year limited warranty

Model Numbers

46070V Solid State Portable Headlight with Direct Power Supply/Charger

Accessory Numbers

05232-U Soft Case

Replacement Parts

49642V Soft Terry Cloth Headband—case of 10

72270 Rechargeable Battery

74170 Direct Power Supply (U.S.)



For international part numbers please contact Welch Allyn Customer Service or your local Welch Allyn distributor.