



Welch Allyn KleenSpec® Disposable Vaginal Specula and Cordless Illumination System

An advancement in women's health technology
that's nothing short of brilliant.



WelchAllyn®

Advancing Frontline Care™

Now you can have better visualization of the exam area without the hassle of cords.

The Welch Allyn KleenSpec® Cordless Illumination System features a cordless design, so you'll never have to worry about cords getting in your way or breaking again.

Plus, with 5,500° K white light-emitting diodes (LEDs) as opposed to halogen, you'll experience a superior whiter, brighter light that provides you with true tissue color during an exam.



LED VS. HALOGEN

The Welch Allyn KleenSpec® Cordless Illumination System will change the way you see your patients and easily fits right into KleenSpec® Disposable Vaginal Speculum.

Cordless and convenient, with you in mind

- No cord to get in the way, clean or break during procedures
- Continuous on time of 80 minutes
- Lithium ion battery provides a recharge time of six hours
- Automatic on/off—insert the illumination system into the vaginal speculum and it automatically illuminates, remove it and the light goes off

Superior light source, for improved patient exams

- No bulb replacement necessary for reduced cost of ownership
- 5,500° K white LED technology for a whiter, brighter light
- Truer tissue rendition and uniform spot—no dark or hot spots—for enhanced visualization of the exam area
- Produces much less heat than traditional lamps for greater comfort



790 Series Illumination System
Shown with #79910 Cordless Illuminator
with Charging Station



The 590 Series KleenSpec® Vaginal Speculum has been redesigned to better meet your needs.

New ergonomic design, to improve patient exams

- A wider handle for improved maneuverability and better balance
- The ability to better manipulate during exams
- A newly designed light pipe that provides a homogeneous spot and eliminates glare back so all of the light is projected forward
- Smooth, molded acrylic that will not pinch or bind tissue for greater patient comfort
- Stronger grip, and the yoke is more flexible than low-cost substitutes

Disposable, for your convenience

- Single-use speculum eliminates the time and expense of reprocessing metal specs
- Improves patient safety by reducing the chance of cross-contamination



#59000

#59001

#59004

590 Series KleenSpec® Vaginal Speculum
Available in small (#59000), medium (#59001) and large (#59004) specula



#78810

The 590 Series KleenSpec® Vaginal Speculum is compatible with the corded illumination system (#78810) as well as the new cordless illumination system (#79910).

Ordering Information

KleenSpec® Vaginal Specula—Premium 590 Series

- 59000 KleenSpec® 590 Series Premium Disposable Vaginal Specula, Small (24/box; 4 boxes/case)
- 59001 KleenSpec® 590 Series Premium Disposable Vaginal Specula, Medium (24/box; 4 boxes/case)
- 59004 KleenSpec® 590 Series Premium Disposable Vaginal Specula, Large (19/box; 4 boxes/case)



KleenSpec® Vaginal Specula—Premium 590 Series with Smoke Tube

- 59005 KleenSpec® 590 Series Premium Disposable Vaginal Specula with Smoke Tube, Small (11/box; 4 boxes/case)
- 59006 KleenSpec® 590 Series Premium Disposable Vaginal Specula with Smoke Tube, Medium (11/box; 4 boxes/case)



KleenSpec® Cordless Illumination System

- 79910 Complete KleenSpec® Cordless Illuminator with Charging Station
- 79915 Same as above, for export purposes—specify country and voltage; conforms to IEC 601-1
- 79920 Complete KleenSpec® Cordless Illuminator with Charging Station and Auxiliary Power Cord
- 79925 Same as above, for export purposes—specify country and voltage; conforms to IEC 601-1
- 79900 KleenSpec® Cordless Illuminator
- 73910 Charging Station for Cordless Illuminator; For export purposes: 73912 (European), 73914 (UK), 73916 (Australia)
- 75910 Auxiliary Power Cord Accessory for Cordless Illuminator; For export purposes: 75912 (European), 75914 (UK), 75916 (Australia)
- 59010 Disposable Sheath for Cordless Illumination System (100/box; 5 boxes/case)



In addition, we offer a complete line of Vaginal Specula and accessories to meet your needs:

KleenSpec® Vaginal Specula—Premium 580 Series with Sheath

- 58000S KleenSpec® 580 Series Premium Disposable Vaginal Specula with Sheath, Small (20/box; 4 boxes/case)
- 58001S KleenSpec® 580 Series Premium Disposable Vaginal Specula with Sheath, Medium (20/box; 4 boxes/case)
- 58004S KleenSpec® 580 Series Premium Disposable Vaginal Specula with Sheath, Large (15/box; 4 boxes/case)



KleenSpec® Vaginal Specula—Economy Series

- 58600 KleenSpec® 586 Series Economy Disposable Vaginal Specula, Small (25/box; 4 boxes/case)
- 58601 KleenSpec® 586 Series Economy Disposable Vaginal Specula, Medium (25/box; 4 boxes/case)
- 78600 Expendable Illuminator (6/box; 6 boxes/case)—Can be sterilized 50 times before disposal

KleenSpec® Corded Illumination System

- 78810 KleenSpec® Illumination System (with cord)—4.7 V Halogen Light Handle with on/off switch and transformer. 110-130 V, 60 Hz for North America. UL Listed, CSA Certified
- 78815 Same as 78810, for export purposes—Specify country and voltage; conforms to IEC 601-1
- 78800 Permanent 4.7 V Halogen Light Handle with on/off Switch
- 73410 Transformer with 5' Cord—110-130 V, 50 Hz. UL Listed, CSA Certified
- 73315 Same as 73410, but for export purposes—Specify country and voltage; conforms to IEC 601-1
- 08800-U 4.6 V to 4.7 V Halogen Replacement Lamp for 78810 Illumination System
- 07800-U 6 V Halogen Replacement Lamp for 78010 Illumination System
- 58010 Disposable Sheath for Corded Illumination System (25/box; 5 boxes/case)



KleenSpec® Storage Systems

- 78820 KleenSpec® Illumination Storage System
- 58014 KleenSpec® Vaginal Speculum Dispenser

Green Series™ Exam Light IV— Another LED innovation from Welch Allyn.

- 48810 GS Exam Light IV with Table/Wall Mount
- 48950 Mobile Stand for GS Exam Light IV
- 48955 Table/Wall Mount for GS Exam Light IV
- 52640 Disposable Sheath for Exam Light III and GS Exam Light IV (25/box, 5 boxes/case)



For more info, go to www.welchallyn.com/greeneries

4341 State Street Road, PO Box 220, Skaneateles Falls, NY 13153-0220 USA
(p) 800.535.6663 (f) 315.685.3361 www.welchallyn.com



WelchAllyn®

Advancing Frontline Care™