



Introducing the new KleenSpec® Disposable Vaginal Specula
and Cordless Illumination System



An advancement in Women's Health
technology that's nothing short of brilliant.

WelchAllyn®

Advancing Frontline Care™

Now you can have better visualisation of the exam area without the hassle of cords.

The new KleenSpec® Disposable Vaginal Speculum with the innovative Cordless Illumination System offers greater convenience for you, added comfort for your patients.

The Welch Allyn KleenSpec® Vaginal Specula and Illumination System has been totally redesigned to provide the key features you've come to expect – plus a number of important and exciting enhancements.

With so many new time-saving features and benefits, the new KleenSpec® Disposable Vaginal Specula and Cordless Illumination System are a must for every office and hospital.




The new 590 Series KleenSpec® Vaginal Speculum has all of the quality features you rely upon and more.

- > **Manoeuvrability** – new wider handle provides improved ergonomics, better balance, and it is much easier to manipulate during exams.
- > **Superior Illumination** – newly designed light pipe provides a homogeneous spot and eliminates glare back – all of the light is projected forward right where you need it.
- > **Comfortable** – smooth, moulded acrylic will not pinch or bind tissue for greater patient comfort during the exam.
- > **Disposable** – Welch Allyn's single-use speculum eliminates the time and expense of reprocessing and reduces the chance of cross-contamination.
- > **Safe and Dependable** – stronger grip and the yoke is more flexible than low-cost substitutes.
- > **Compatibility** – the new 590 Series KleenSpec® Vaginal Speculum is compatible with the current illumination system (78812/78814) as well as the new illumination system (79912/79914).



590 Series KleenSpec® Vaginal Speculum
Available in small (59000), medium (59001) and large (59004) specula



The new Welch Allyn KleenSpec® Disposable Vaginal Specula 790 Series Illumination System will change the way you see your patients.

- > **Superior Visualisation** – advanced LED light source provides better visualisation of the exam area by supplying brighter, whiter light, right where it is needed.
- > **Cordless and Rechargeable** – no cord to get in the way or break during procedures, continuous on time of 80 minutes, recharge time of six hours, lithium ion battery technology. The new cordless illumination system easily fits right into the new KleenSpec® Disposable Vaginal Speculum.
- > **Automatic On/Off** – no more hassle with turning the unit on and off. When the new KleenSpec® Illumination System is inserted into the 590 Vaginal Speculum it is automatically illuminated, and once it is removed the light goes off.
- > **Versatility** – use your new Illumination System as an all-purpose light source as well as with the KleenSpec® Disposable Vaginal Speculum.
- > **Convenient** – no cords to clean, much less heat than traditional lamps.



790 Series Illumination System

Shown with 79922/79924 KleenSpec® Disposable Vaginal Specula Cordless Illuminator with Charging Station and Auxiliary Power Cord

The new KleenSpec® Disposable Vaginal Specula and Cordless Illumination System

Ordering Information

KleenSpec® Vaginal Specula – Premium 590 Series

- 59000 KleenSpec® 590 Series Premium Disposable Vaginal Specula, Small (24/box; 4 boxes/case)
- 59001 KleenSpec® 590 Series Premium Disposable Vaginal Specula, Medium (24/box; 4 boxes/case)
- 59004 KleenSpec® 590 Series Premium Disposable Vaginal Specula, Large (19/box; 4 boxes/case)

KleenSpec® Box Series

- 59000-B KleenSpec® 590 Series Premium Disposable Vaginal Specula, Small (24/box)
- 59001-B KleenSpec® 590 Series Premium Disposable Vaginal Specula, Medium (24/box)
- 59004-B KleenSpec® 590 Series Premium Disposable Vaginal Specula, Large (19/box)

KleenSpec® Vaginal Specula – Premium 590 Series w/Smoke Tube

- 59005 KleenSpec® 590 Series Premium Disposable Vaginal Specula with Smoke Tube, Small (12/box; 4 boxes/case)
- 59006 KleenSpec® 590 Series Premium Disposable Vaginal Specula with Smoke Tube, Medium (12/box; 4 boxes/case)

KleenSpec® Disposable Vaginal Specula Illumination System

- 79912 Complete KleenSpec® Disposable Vaginal Specula Cordless Illuminator with Charging Station (Europe)
- 79914 Complete KleenSpec® Disposable Vaginal Specula Cordless Illuminator with Charging Station (UK)
- 79922 Complete KleenSpec® Disposable Vaginal Specula Cordless Illuminator with Charging Station and Auxiliary Power Cord, conforms to IEC 601-1 (Europe)
- 79924 Complete KleenSpec® Disposable Vaginal Specula Cordless Illuminator with Charging Station and Auxiliary Power Cord, conforms to IEC 601-1 (UK)
- 79900 KleenSpec® Disposable Vaginal Specula Cordless Illuminator
- 73912 Charging Station for Cordless Illuminator (Europe)
- 73914 Charging Station for Cordless Illuminator (UK)
- 75912 Auxiliary Power Cord Accessory for Cordless Illuminator (Europe)
- 75914 Auxiliary Power Cord Accessory for Cordless Illuminator (UK)



In addition, we offer a complete line of Vaginal Specula and Accessories to meet your needs:

KleenSpec® Vaginal Specula – Premium 580 Series w/Sheath

- 58000S KleenSpec® 580 Series Premium Disposable Vaginal Specula with Sheath, Small (20/box; 4 boxes/case)
- 58001S KleenSpec® 580 Series Premium Disposable Vaginal Specula with Sheath, Medium (20/box; 4 boxes/case)
- 58004S KleenSpec® 580 Series Premium Disposable Vaginal Specula with Sheath, Large (15/box; 4 boxes/case)

KleenSpec® Box Series

- 58000S-B KleenSpec® 580 Series Premium Disposable Vaginal Specula w/Sheath, Small (20/box)
- 58001S-B KleenSpec® 580 Series Premium Disposable Vaginal Specula w/Sheath, Medium (20/box)

KleenSpec® Vaginal Specula – Economy Series

- 58600 KleenSpec® 586 Series Economy Disposable Vaginal Specula, Small (25/box; 4 boxes/case)
- 58601 KleenSpec® 586 Series Economy Disposable Vaginal Specula, Medium (25/box; 4 boxes/case)
- 78600 Expendable Illuminator (6/box; 6 boxes/case) – Can be sterilized 50 times before disposal

KleenSpec® Disposable Vaginal Specula Illumination System

- 78812 KleenSpec® Illumination System – 4.7 V Halogen Light Handle with on/off switch and transformer (Europe)
- 78814 KleenSpec® Illumination System – 4.7 V Halogen Light Handle with on/off switch and transformer (UK)
- 78800 Permanent 4.7 V Halogen Light Handle with On/Off Switch
- 73412 Transformer with 1.5m Cord – UL Listed, CSA Certified, conforms to IEC 601-1 (Europe)
- 73314 Transformer with 1.5m Cord – UL Listed, CSA Certified, conforms to IEC 601-1 (UK)
- 08800-U 4.6 V to 4.7 V Halogen Replacement Lamp for 78812/78814 Illumination System
- 07800-U 6 V Halogen Replacement Lamp for 78012/78014 Illumination System

KleenSpec® Storage Systems

- 78820 KleenSpec® Illumination Storage System
- 58014 KleenSpec® Vaginal Speculum Dispenser

Patents Pending



UK Customer Services

Welch Allyn UK Ltd, Cublington Road,
Aston Abbots, Buckinghamshire HP22 4ND, England
Tel: +44 (0) 207 365 6780
Fax: +44 (0) 207 365 9694

Welch Allyn South Africa (Pty) Ltd, 394 Surrey Avenue,
Ferndale, Randburg, 2194, South Africa
P O Box 809, Fountainsbleau, South Africa 2032
Tel: +27 (0) 11 777-7555
Fax: +27 (0) 11 777-7556

Welch Allyn Customer Services
Navan Business Park, Dublin Road,
Navan, Co. Meath, Republic of Ireland
Tel: +353 46 90 67790
Fax: +353 46 90 67754

WelchAllyn®

Advancing Frontline Care™