

**Parts List For Welch Allyn AED20 and PIC Defibrillators
Available from ZOLL Medical – April 2009**

Product	Part#	Additional Versions	Description
PIC	001515-U		Electrode: Internal Electrodes: 25mm (1") - (Set of 2)
PIC	001516-U		Electrode: Internal Electrodes: 50mm (2") - (Set of 2)
PIC	001537-U		Electrode: Pediatric External Paddle Adapter - For Use With 971108 External Defibrillation Paddles (Set of 2)
PIC	001647-U		Power: Welch Allyn SuperPac Nickel Metal Hydride (NiMH) Rechargeable Battery
PIC	001739-U		Paper: ECG Plain White Paper - 50mm - (3 rolls =one box)
PIC	001790-U		Cable: PIC 3-Lead ECG (AHA)
PIC	001794-U		Cable: 3LD PATIENT CBLE INTL PIC
PIC	001795-U		Cable: PIC 5-Lead ECG (AHA)
PIC	001796-U		Cable: 5LD PATIENT CBLE INTL PIC
PIC	001798-U		Cable: PIC 3-lead ECG (AHA) for use with 12-lead PIC50
AED20	001829-E	& 6 other languages -C,-D,-F,-G,-I,-S	Power: AED 20 Rechargeable (NiMH) Battery
AED20	00183-0E	& 6 other languages as 001830- -C,-D,-F,-G,-I,-S	Power: AED 20 Non-Rechargeable (Lithium) Battery
PIC	001835-U		Cable: Welch Allyn "Breakaway" (AHA) Trunk Cable and 4 leadwire set only
PIC	001836-U		Cable: Detachable 6 lead-wire set (AHA) for 001834 "Breakaway" 12-Lead (AHA) Patient Cable – 6 lead pre-cordial set only
PIC	001838-U		Cable: PTNT CABLE 4LD 12POLE IEC
PIC	001839-U		Cable: PATIENT CABLE V LEADS IEC
PIC/AED20	00185-3		Electrode: "Leads-In" Defibrillation Pads (1 pouch / 2 pads)
PIC/AED20	001855-U		Electrode: "Leads-In" Defibrillation Pads (1 box / 10 pair)
PIC/AED20	001857-E		Electrode: "Leads-Out" (Pre-Attached) Defibrillation Pads (1 pouch / 2 pads)
PIC/AED20	001858-E		Electrode: Leads-Out (Pre-Attached) Defibrillation Electrodes (1 box / 10 pair)
PIC	001932-U		Cable: 10 foot (3 meter) sensor cable for all disposable sensors
PIC	001934-U		Sensor: Welch Allyn PIC CO2 (Capnography) sensor
PIC	001951-U		Adapter: Adult Disposable Airway Adapters - (1 box / 50 adapters)
PIC	001954-U		Cable: Welch Allyn PIC To CO2 Sensor Adapter (pink connector)
PIC	001955-U		Cable: Welch Allyn PIC to Nellcor Brand Sensor Adapter Cable (green connector)
PIC	001956-U		Cable: Welch Allyn PIC temperature adapter cable - sensor / probe cable (001932) and sensors acquired separately
PIC	001957-U		Cable: Invasive Blood Pressure Adapter - Blue (Welch Allyn PIC IBP Connection Point to Interface / Transducer Cable)
AED20	001962-U		Reporting: SmartLink Event Pro Event Reporting Software (for use with AED20 and AED10)
PIC	001966-E		Cable: Welch Allyn PIC to Nellcor Brand Sensor Adapter (gray connector)
PIC	001971-E		Adapter: PIC Serial Cable Bracket/Gender Changer (for use with SmartLink Wireless Data Transition Software)
PIC	001974-U		Hose: Welch Allyn PIC - Non-Invasive Blood Pressure Cuff Hose - 10 foot (3 meter) extension hose for use with Welch Allyn Single Hose Blood Pressure Cuffs with Female Locking Connector (5082-xx-3MQ style cuffs)*.
AED20	002120-U		Cable: AED 20 Data Transmission Cable for use with SmartLink Event Pro or Lite Reporting Software Program (AED20 Only)

**Parts List For Welch Allyn AED20 and PIC Defibrillators
Available from ZOLL Medical On April 27, 2009**

Product	Part#	Additional Versions	Description
AED20	002128-U		Cable: AED 20 - Single Channel Monitoring ECG (Lead II – Using IEC Identifiers)
AED20	00213-0		Cable: AED 20 - Single Channel Monitoring ECG (Lead II – Using AHA Identifiers)
AED20	002174-U		Adapter: Pediatric Energy Reducer - for use with AED 20 only
AED20	0-02181		Mount: AED 20 White L-shaped Metal Wall Bracket -- includes mounting hardware
PIC	560640		Lens: Welch Allyn PIC Cover (CLEAR) - for rubberized front panel PIC units
PIC	570643		Power: Invasive Blood Pressure Auxiliary Power Adapter - Required for IBP Feature Function (Part of 971017 kit)
PIC	570663		Case: PIC Ambulance Mount Carrying Case - Required For Use With Ambulance Mount (981114*) This PIC carrying case is specifically designed for use with the Ambulance Mount (included integral charging ports for direct connection to Ambulance Mount)
PIC	90021-4		Case: PIC Carrying Case (includes: top pouch, side pouches, cable management straps)
PIC	971029-E	& 9 other languages A, F, G, I, M, S, T,U, W	Power: Paddle Tray / Integral Charger / Bustle Pack (using AC Power) for PIC - paddles purchased separately
PIC	971107		Cable: PIC Hands Free Defibrillation / Pacing Cable
PIC	971108-E	& 5 other languages F, G, I, PL and S	Electrode: Deluxe External Paddles With Energy Select, Charge Buttons and Energy Discharge Buttons On Paddles - includes one tube of defibrillation gel
PIC	981114		Mount: Ambulance Mount - Secures PIC to vehicle and trickle charges battery using 12V power source
PIC	981117E		Power: PIC Basic Power Supply (110V/220V) - Power supply to operate PIC and charge SuperPac battery while in the device.
PIC	981118		Power: PIC Ambulance Mount Power Supply (110V/220V) - Power supply to operate PIC and charge SuperPac battery while in the device.
PIC	981122E	& 4 other languages A,T,U,W	Power: Dual Slot PIC SuperPac Battery Charger / includes battery reconditioned (charges 2 PIC SuperPac Batteries)
PIC	981125E	& 4 other languages A,T,U,W	Power: Single Slot AED 20 Battery Charger
AED20	810-2484-02		Manual: AED20 Multilingual User Manual and In-service CD/DVD
PIC	810-2558-01		Manual: PIC Multilingual User Manual CD (PIC30/PIC40/PIC50)

Language / Mains Power Destination Codes:

E-English / US, **A**-English / Australia, **U**-English / UK, **T**-English / SoAfrica, **W**-English / Europe, **S**-Spanish / Spain, **G**-German / Germany, **F**-French / France, **I**-Italian / Italy, **C**-Mandarin Chinese / China, **D**-Danish / Denmark, **PL**-Polish, **M**-Mexico

To Place An Order, Contact Your Authorized ZOLL Dealer or ZOLL Medical Corporation