


Welch Allyn RetinaVue™ Network


Diabetic retinal exams made **SIMPLE** and **AFFORDABLE** for primary care



Diabetic retinopathy is the leading cause of blindness among working-age adults.



Documented compliance with annual retinal exams is often **less than 50%**.



With early detection, **95% of vision loss cases can be prevented.**



Now you can easily detect diabetic retinopathy in primary-care settings using a **SIMPLE and EFFECTIVE** four-step process.



Welch Allyn RetinaVue™ 100 Imager

BREAKTHROUGH CAMERA TECHNOLOGY IS LOW COST AND EASY TO USE

- One-third the cost of desktop fundus cameras
- Touchless auto-focus and image capture support ease of use
- Automatic evaluation of image helps ensure quality
- Fast, non-mydratic exam promotes patient comfort

RetinaVue Network is designed for healthcare providers in primary care settings who want to improve management of their patients with diabetes.

Acquire

Acquire non-mydratiatic fundus images **in minutes** with the simple, automated Welch Allyn RetinaVue 100 Imager

Transfer

Encrypted images sent via secure **HIPAA-compliant** RetinaVue Network Software

Analyze

Images are evaluated by a **board-certified ophthalmologist**

Report

Diagnostic report and referral/screening plan **in one business day**
—usually in 90 minutes


PATIENT COMPLIANCE RATES CAN DOUBLE IN 12 MONTHS

Single Primary Care Clinic Retinal Assessment Compliance Results using the RetinaVue Network

Garg S, Jani PD, Kshirsagar AV, King B, Chaum E. Telemedicine and Retinal Images for Improving Diabetic Retinopathy Evaluation. Arch Intern Med. 2012; Oct 1: 1-2.

* Patients with a Documented Retinal Assessment, %


COMPLIANCE STUDY CHART



MOST TOP HEALTHCARE PLANS REIMBURSE FOR SCREENINGS

Four of the top five commercial healthcare plans provide coverage for RetinaVue diabetic retinopathy screening. Consult your payer partners to understand coding and coverage options.

¹ CPT ©2015, American Medical Association



IMPROVE QUALITY METRICS AND REVENUE

Annual retinal examinations are included in the NCQA HEDIS ratings program, the Medicare Advantage STAR quality rating program and Medicare quality rating programs.





Diagnostic interpretation services are a required component of the RetinaVue Network and are contracted directly with RetinaVue P.C.



Diagnostic interpretation and electronic report turned around in one business day. Average time: 90 minutes

Diagnostic report includes retinal images and contains all relevant ICD codes

Referral/Screening plan clearly details next steps

Comments, signature and license number from board-certified ophthalmologist



RetinaVue™ Network
Diagnostic Report

Patient Information

Name: Smith, Mary
Medical Record Number: 512124
Date of Birth: 8/14/1970
Referring Physician: William Brown, MD
Referring Clinic: Primary Care Associates

Retinal Image Assessment and Management Plan

Fundus Photograph of Left Eye (OS):



Diagnosis for Left Eye (OS):
NPDR Moderate with CSME

Fundus Photograph of Right Eye (OD):



Diagnosis for Right Eye (OD):
NPDR Severe with CSME

ICD-10 Diagnosis Codes:

E11.331 Type 2 diabetes mellitus with moderate nonproliferative diabetic retinopathy with macular edema
E11.341 Type 2 diabetes mellitus with severe nonproliferative diabetic retinopathy with macular edema.

Recommendation and Management Plan:

Refer to a retina specialist for evaluation of evidence of severe diabetic eye disease

Comments:
None

Reviewing Physician: John Doe, MD
TN License Number: 654321
Report Generated on: 12/16/2015

Diagnosis and Management Plan
RetinaVue Network - Signature Panel

This report is provided by RetinaVue, P.C. through its licensed and board-certified ophthalmologists. This report specifies the presence of disease in the posterior pole of the retina for diagnostic, management and referral purposes. It does not take place of a regular eye examination for the purpose of assessing the presence of glaucoma, cataract, anterior segment diseases, peripheral retina disease or other possible vision threatening conditions. For customer support, please call 1.800.535.6663, 9 a.m. to 5 p.m. EST.
© 2015 Welch Allyn MC13062

Ordering Information



RV100-B*
RetinaVue 100 Imager



723427*
Topcon® TRC NW-400
Retinal Camera



723480
K2 Camera Table



723481
Lenovo® Thinkpad

* RetinaVue Network software is shipped with either camera. System requirements: Windows 7, Windows 8.1, Windows 10 on latest SP; RAM 2 GB min; CPU 2 GHz or equivalent; 150 MB free HDD space; high speed internet connection; USB port (for RV100-B); Ethernet port (for 723427)

Welch Allyn **VISION FOR ALL**

Delivering simple and affordable screening solutions in primary care settings to help eradicate some of the leading causes of preventable blindness.

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361

