

WelchAllyn®



Welch Allyn  
Connex® Spot Monitor

# Welch Allyn Connex Spot Monitor

Accurate, automatic vitals sent to your EMR from the point of care to help enhance staff productivity, patient outcomes and risk reduction

Choose the Connex Spot Monitor for wireless data transfer, accurate vital signs measurement across your patient population and compatibility with your EMR—now and in the future.

## Accurate readings across your patient population

- > Help ensure accurate hypertension diagnosis without slowing down your workflow with unique SureBP® technology and blood pressure averaging.
  - Research shows that more than half of patients could be misdiagnosed in a hypertensive category using the single-reading method.<sup>1</sup>
- > Provide SpO<sub>2</sub> accuracy in challenging conditions including motion, thick fingers, dark skin tone and low perfusion with Masimo® SET® and Nonin® PureSAT® pulse oximetry technology.<sup>2</sup>

## Improve data accuracy, efficiency and mobility

- > The Connex Spot Monitor is designed to help eliminate data entry errors in the patient chart with wireless transfer of vital signs.
- > Increase staff efficiency by eliminating typing into patient records.
- > Provide freedom from wires and cables, giving you the flexibility to easily work in limited space.



# Configurable design

allows you to customize the monitor to the needs of your practice and your patients

Every Connex Spot Monitor features an easy-to-use touchscreen display and reliable SureBP technology. You can also select from a number of options to get exactly the monitor you need, including:



## Thermometry

- SureTemp® Plus Thermometry
- Braun ThermoScan® PRO 6000 Ear Thermometry

## Pulse Oximetry

- Masimo SET SpO<sub>2</sub>
- Nonin PureSAT SpO<sub>2</sub>
- Nellcor® OxiMax® SpO<sub>2</sub>

## Mounting and Mobility Options

- **Integrated Wall System** combines physical assessment and digital vital signs for a complete patient exam workflow
- **Mobile Work Surface Stand** provides added work space for note taking, storage or laptop placement
- **Accessory Power Management Stand** features an integrated battery to ensure the device gets you through a busy day
- **Standard Basket Stand** provides the simple portability you need with a minimal footprint



## Connectivity for Ambulatory Care

- > Import vitals data to leading EMR providers—including Allscripts®, athenahealth®, Cerner®, eClinicalWorks®, Epic®, GE® Centricity™, Greenway® (Intergy®), NextGen® and many others—to ensure reliable device connectivity and to help you leverage your EMR investment.
- > Connect your Connex Spot Monitor to scales from Detecto®, Health o meter® and SECA to accurately document height, weight and BMI.

## Protect your investment with service programs

- > **Partners in Care™ Comprehensive Partnership Programs** help minimize the cost and complexity of servicing and maintaining your monitors.
- > **Enterprise Deployment Option** enables your monitors to arrive configured to your specification right out of the box.

<sup>1</sup>Robert Smith, MD. [Blood pressure averaging methodology: Decreasing the rate of misdiagnosing hypertension](#). 2014.

<sup>2</sup>Nonin Medical, [Patient conditions can be unpredictable. Pulse oximeter accuracy shouldn't be](#). 2014.

## Configurations Available

### Connex Spot Monitor

	7100 Series	7300 Series
<b>Parameters</b>		
SureBP technology with BP averaging protocol	●	●
Interval monitoring	●	●
SureTemp Plus thermometry	●	●
Braun ear thermometry	●	●
Masimo pulse oximetry		●
Nellcor pulse oximetry		●
Nonin pulse oximetry	●	●
<b>Communications</b>		
USB	●	●
Ethernet	●	●
Bluetooth®		●

	7100 Series	7300 Series
<b>Accessories</b>		
Classic Basket Stand	●	●
Work Surface Stand	●	●
Advanced Power Management Stand	●	●
GCX Wall Mount	●	●
777 Wallboard Series	●	●
<b>Scale Connectivity</b>		
SECA	USB and Wireless	USB and Wireless
Health o meter®	USB and Wireless	USB and Wireless
Detecto®	USB	USB
<b>Warranty</b>		
Two-year out of box	●	●
Partners in Care services: one, two or five years (optional)	●	●

## ACCESSORIES

<b>7000-APM</b>	Connex Spot Monitor Accessory Power Management Stand
<b>7000-MS3</b>	Connex Spot Monitor Classic Mobile Stand
<b>7000-GCX</b>	Connex Spot Monitor Wall Channel
<b>7000-MWS</b>	Connex Spot Monitor Mobile Work Surface Stand
<b>77794-MCSM</b>	Green Series 777 Integrated Wall System with Coaxial Ophthalmoscope
<b>77794-2MCSM</b>	Green Series 777 Integrated Wall System with PanOptic Ophthalmoscope
<b>7000-PS</b>	Power Supply (required if not ordering a stand or wall mount)
<b>BATT22</b>	Connex Spot Monitor Lithium-Ion Battery
<b>BATT99</b>	Accessory Power Management Stand Lithium-Ion Battery
<b>106275</b>	Connectivity Accessory Kit USB Cable
<b>4500-35</b>	Blood Pressure Hose with FlexiPort® (10 ft)
<b>7000-33</b>	Neonatal Blood Pressure Hose (10 ft)
<b>02895-100</b>	SureTemp Plus Thermometer Rectal Probe and Well
<b>718584</b>	Braun PRO 6000 Ear Thermometer Tether
<b>3278-010</b>	Nonin Reusable Adult Sensor (2 meters)
<b>2360-010</b>	Nonin Reusable Pediatric Sensor (2 meters)
<b>S1-CSM-5</b>	Connex Spot Monitor Comprehensive Partnership Program, 5 years*
<b>S2-CSM-5</b>	Connex Spot Monitor Biomed Partnership Program, 5 years*

\*One- and two-year service plans also available.



Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA  
(p) 800.535.6663 (f) 315.685.3361

[WWW.WELCHALLYN.COM](http://WWW.WELCHALLYN.COM)



**WelchAllyn®**