



Welch Allyn  
CP 150™ Electrocardiograph  
with Optional Spirometry  
Simple. Fast. Connected.



Optional wireless connectivity and DICOM interoperability

# Help your practice run more efficiently with the Welch Allyn CP 150 Resting ECG—now with optional wireless capability.

Simple, streamlined ECG testing designed to help you improve workflows and get more from your investments


The CP 150 combines the best traditional box ECG features (e.g., full-size printouts and clinical decision support with MEANS interpretation) with connectivity enhancements designed especially for today's primary care practices.

In fact, **94% of CP 150 users surveyed** say the touchscreen display and simple user interface have helped them improve their practice efficiency.\*

\* Source: Survey by Walker Information, Inc. of clinicians who have purchased the CP 150.


94%

## Simple



### Intuitive touchscreen display

Can improve accuracy, minimize training and save time.



### Lead quality display


Help avoid repeat ECGs by identifying lead connection problems before completing a test.

## Fast




### Instant-On Feature

Help eliminate typical start-up delays.



### 12-Lead ECG Preview

Large display provides three, six or 12-lead ECG preview options.



### ECG Print Preview

Can save time and reduce paper costs by verifying the ECG test quality prior to printout.



Send ECG results to:



Network folder as PDF  
(For import into EMR)



DICOM® image server  
(EMR, CVIS, PACS)

## Connected

### Wireless connectivity increases device mobility and utilization

- Wireless connectivity through an integrated 802.11 radio
- Share results immediately to put the data where you need it—in your EMR/CVIS/PACS—right away
- Meets globally accepted data encryption standards for enhanced patient information security



### Connectivity options tailored to enhance your practice's workflow efficiency

#### Wireless PDF transfer to network folders

#### DICOM® conformance enhances interoperability with multiple DICOM-compatible platforms

- Bidirectional workflow with your EMR/CVIS<sup>1</sup> or PACS<sup>2</sup> without the middleware
- Reduce patient data re-entry by importing the DICOM modality worklist
- Upgrade your existing CP 150 to make it DICOM compatible

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

<sup>1</sup> Cardiovascular information system

<sup>2</sup> Picture archiving and communication system

# Spirometry Option

The CP 150 lets you perform two reimbursable tests with one device.

- Perform FVC and pre- and post-BD testing with automatic interpretation and reversibility analysis
- Customize report formats and screen displays to fit your needs and workflow
- Integrated ambient sensors can eliminate the need for the user to find and enter values for temperature, pressure and humidity
- Manage patients with COPD to help you maximize HEDIS scores



The CP 150 spirometry option features four incentive screens to help educate and coach your patients:



## Ordering Information

### Selected Device Configurations\*:

CP150-IENB	CP 150 Non-Interpretive Electrocardiograph
CP150A-IENB	CP 150 Interpretive Electrocardiograph
CP150AW-IENB	CP 150 Interpretive Electrocardiograph with Wireless Capability
CP150AD-IENB	CP 150 Interpretive Electrocardiograph with DICOM
CP150AS-IENB	CP 150 Interpretive Electrocardiograph with Spirometry
CP150AWD-IENB	CP 150 Interpretive Electrocardiograph with Wireless Capability and DICOM
CP150AWDS-IENB	CP 150 Interpretive Electrocardiograph with Wireless Capability, DICOM and Spirometry

### Accessories:

715006	ECG Multifunction Electrode Adaptor
45008-0000	Resting Tab Electrodes (box of 1000)
105353	CP 150 Recording Paper (case, 1000 sheets)
719653	10-lead Patient Cable, AHA, Banana, CP 150
105341	CP 150 Office Cart
105342	CP 150 Hospital Cart
105343	CP 150 Cable Arm & Shelf Cart

\*For more device configurations, see CP 150 Part Number Configuration document.

### Spirometry Accessories:

105660	Spirometry Upgrade for CP 150 (works with any standard CP150A-IENB or CP150-IENB configuration)
720705	Disposable Flow Transducers, 25 pk, CP 150
720706	Disposable Flow Transducers, 100 pk, CP 150 / Spiro

### Partners in Care<sup>SM</sup> Service Options:

S1-CP150	CP 150 Comprehensive Partner Program
S1-CP150-2	CP 150 Comprehensive Partner Program, 2 YR
S1-CP150S	CP 150 Spiro Comprehensive Partner Program
S1-CP150S-2	CP 150 Spiro Comprehensive Partner Program, 2 YR

### Connectivity Upgrades:

406814	CP 50/150 Connectivity Kit
106736	CP 150 DICOM Upgrade Kit

### Warranty Information:

CP 150 Electrocardiograph	—3 years
Patient Cable/Battery	—90 days
Spirometer	—1 year

Contact your Welch Allyn representative or visit [www.welchallyn.com](http://www.welchallyn.com) for more information.

Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA  
(p) 800.535.6663 (f) 315.685.3361

[WWW.WELCHALLYN.COM](http://WWW.WELCHALLYN.COM)

© 2016 Welch Allyn MC13436 80021032 Ver A 2016-11

**WelchAllyn**<sup>®</sup>