





Multi-parameter Vital Signs Devices



	 Spot Vital Signs® Device	 Connex® Spot Monitor (CSM)	 Connex® Vital Signs Monitor (CVSM)	 Connex® Integrated Wall System (CIWS)
Workflows	Spot	Spot, Office, Intervals	Spot, Office, Intervals, Continuous	Spot, Office, Intervals, Continuous
Device Type	Basic	Intermediate	Advanced	Advanced
Patient Populations	Adult, Pediatric	Adult, Pediatric, Neonatal	Adult, Pediatric, Neonatal	Adult, Pediatric, Neonatal
Device Alarms		Yes	Yes	Yes
Third-Party Alarm Notification			Yes	Yes
Parameters	Blood Pressure, Temperature, SpO ₂	Blood Pressure, Temperature, SpO ₂	Blood Pressure, Temperature, SpO ₂ , etCO ₂ , SpHb, ECG, RRa	Blood Pressure, Temperature, SpO ₂ , SpHb, ECG
SureBP® Technology (15-sec)		Yes	Yes	Yes
Blood Pressure Averaging		Yes	Yes	Yes
Contact-free Monitoring			Respiration, Bed Exits, Turn Reminders	
Diagnostic Handles				Otoscope, Ophthalmoscope
SpO₂ Options	Nellcor®, Masimo®	Nellcor, Masimo, Nonin®	Nellcor, Masimo	Nellcor, Masimo
Temperature Options	SureTemp® Technology	SureTemp Plus, Braun PRO 6000	SureTemp Plus, Braun PRO 6000	SureTemp Plus, Braun PRO 6000
Early Warning Scoring		Up to 6 protocols	1 protocol	1 protocol
Custom Documentation Fields		Up to 20	Up to 20	Up to 20
Wireless Connectivity		WiFi®, Bluetooth®	WiFi	WiFi
Wired Connectivity		USB, Ethernet	USB, Ethernet	USB, Ethernet
EMR Client Support		Thin, Thick	Thin, Thick	Thin, Thick
Display	LED	7" Color LCD, Touchscreen	8.9" Color LCD, Touchscreen	8.9" Color LCD, Touchscreen
Scale Integration		Yes	Yes	Yes
Mounting Options	Basic Stand Wall Mount	Power Management, Stand, Classic Stand, Work Surface Stand, Wall Mount, Integrated Wall System, Desktop Stand	Cable Management, Stand, Basic Stand, Wall Channel	Wall Mounted Unit
Accessory Options		Scales, Barcode, RFID	Scales, Barcode	Scales, Barcode
Memory (# Readings)	1	400	400	400
Launch Year	2001	2015	2010	2011