

# Welch Allyn Universal ECG Cart for Outpatient Care Specifications

When sharing an ECG device across multiple exam rooms, a mobile cart enhances clinical workflow efficiency.

## Universal ECG Cart for Outpatient Care

Part number 9911-024-52

- 1 Work surface is 22" (55.58 cm) x 16" (40.64 cm)
- 2 Two fixed and two swivel casters with brakes
- 3 Accessory bin and cable management
- 4 Cable holster for WAM™ Wireless Acquisition Module or AM12™ patient cable

Depth: 30"  
(76.2 cm)



The Universal ECG Cart for Outpatient Care is designed for use with the following resting ECGs:



ELI™ 280 Resting ECG



ELI 150c Resting ECG



ELI 250c Resting ECG

Contact your representative or visit [www.welchallyn.com](http://www.welchallyn.com) to learn more.



**Clinical excellence. Connected solutions. Continuous innovation.**  
Welch Allyn Cardiology is proud to be powered by Mortara.