



Welch Allyn
SureTemp® Plus Electronic Thermometers



Flexible and accurate technology featuring oral, axillary and rectal modes to meet all your needs.

Enhance a routine task with quick and accurate predictive temperature readings **you can rely on.**



VERSATILE & FAST

- Oral, rectal and axillary modes, to measure all patient types
 - Oral readings: approximately 4- to 6-seconds
 - Pediatric Axillary and Rectal readings: approximately 10- to 13-seconds
 - Adult Axillary readings: approximately 12- to 15-seconds
- Dedicated button helps you quickly switch between modes
- C/F button easily toggles between Celsius and Fahrenheit
- Monitor mode option helps to measure thermal steady state



EASY & INTUITIVE

- Large LCD display features flashing icons
- Recall button displays last temperature taken
- SureTemp Plus 692 model features a 60-second timer to support pulse and respiration readings

SAFE & SECURE

- Soft, single-use disposable probe covers help minimize cross-contamination risk between patients
- Interchangeable probe wells are easy to remove and clean as you switch between modes
- SureTemp Plus 692 model designed with two electronic security features to help protect against theft

Welch Allyn SureTemp® Plus 690 Electronic Thermometer

IDEAL CARE SETTINGS

- Family Practice
- Pediatrics
- Long-term Care
- Internal Medicine
- OB/Gyn



The SureTemp Plus 690 with optional wall mount has storage space for 50 probe covers

Removable probe well for easy cleaning

Large, easy-to-read LCD

C/F button for toggling between Fahrenheit and Celsius

Recall button displays last temperature taken

Convenient storage for 25 probe covers

Battery life indicator

Mode button

Rubberized grips



Welch Allyn SureTemp® Plus 692 Electronic Thermometer

IDEAL CARE SETTINGS

- Urgent Care
- MedSurg
- Maternity
- Emergency Medicine
- Oncology
- Cardiology



The SureTemp Plus 692 with standard wall mount has storage space for 50 probe covers.

Interchangeable rectal probe and well reduce risk of cross-contamination



User-selectable icon and/or words for identification of thermometer mode



Battery life indicator

Removable probe well for easy cleaning

Large, easy-to-read LCD

NEW Backlit display helps in low-light environments

C/F button for toggling between Fahrenheit and Celsius

Recall button displays last temperature taken

Rubberized grips

Convenient storage for 25 probe covers

Pulse timer button/monitor mode

NEW feature makes it easier to switch to Monitor Mode for measuring a patient's thermal steady state (approx. 3 to 5 minutes)

Mode button for selecting oral, pediatric or adult axillary and rectal modes



SureTemp Plus 690, 692 Electronic Thermometer Specifications

Calibration Accuracy	±0.2 °F (±0.1 °C) (Monitor Mode)	
Temperature Predict Time	Oral: Adult Axillary (18 years and older): Pediatric Axillary (17 years and younger): Rectal:	Approximately 4-6 seconds Approximately 12-15 seconds Approximately 10-13 seconds Approximately 10-13 seconds
Unit Dimensions	8.46" x 3.18" x 2.43" (215 mm x 81 mm x 62 mm)	
Weight	12 oz (357 g) (model 690)	12.6 oz (357 g) (model 692)
Ambient Temperature Range	50.0 °F (10 °C) to 104 °F (40.0 °C)	
Patient Temperature Range	80.0 °F (26.7 °C) to 110 °F (43.3 °C)	
Power Source	(3) 1.5 Vdc AA batteries	
Battery Life	Approximately 6,000 temperatures (model 690)	Approximately 5,000 temperatures (model 692)

Ordering Information:

690 MODELS

- 01690-200** SureTemp Plus 690 Electronic Thermometer, Oral Probe and Oral Probe Well with 4ft cord
- 01690-201** SureTemp Plus 690 Electronic Thermometer, Rectal Probe with Rectal Probe Well with 4ft cord
- 01690-300** SureTemp Plus 690 Electronic Thermometer, One per room, Wall Mount, Oral Probe with Oral Probe Well with 9ft cord (must be locked to the wall)
- 01690-301** SureTemp Plus 690 Electronic Thermometer, One Per Room, Wall Mount, Rectal Probe with Rectal Probe Well with 9ft cord (must be locked to the wall)
- 01690-400** SureTemp Plus 690 Electronic Thermometer, Wall Mount, Oral Probe with Oral Probe Well with 4ft cord
- 01690-700** SureTemp Plus 690 Electronic Thermometer, Wall Mount, Oral Probe with Oral Probe Well with 9ft cord with Rolling Stand and Holder for Spare Probe & Well

692 MODELS

- 01692-200** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 4ft cord
- 01692-201** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Rectal Probe and Rectal Probe Well with 4ft cord
- 01692-300** SureTemp Plus 692 Electronic Thermometer, One per room, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord (must be locked to the wall)
- 01692-301** SureTemp Plus 692 Electronic Thermometer, One Per Room, Wall Mount, Rectal Probe and Rectal Probe Well with 9ft cord (must be locked to the wall)
- 01692-700** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord with Rolling Stand and Holder for Spare Probe & Well

ACCESSORIES

- | | |
|--|---|
| 05031-750 Case of 7,500 Disposable Probe Covers | 02893-000 Probe and Well Kit, 4ft, Oral |
| 21326-0000 Wall Holder | 02893-100 Probe and Well Kit, 9ft, Oral |
| 21330-0000 Holder for Spare Probe and Well | 02891-0000 Probe Well, Blue |
| 06138-000 Calibration Key Assembly, 690/692 | 02891-1000 Probe Well, Red |
| 02892-000 Probe and Well Kit, 4ft, Rectal | 08273-000 Thermometer Stand with Holder for Spare Probe and Well |
| 02892-100 Probe and Well Kit, 9ft, Rectal | 21333-0000 Bracket, M600 Stand Upgrade, with Holder for Spare Probe and Well |

Contact your Welch Allyn representative or visit www.welchallyn.com to learn more.



Welch Allyn, Inc.
 4341 State Street Road
 Skaneateles Falls, NY 13153 USA
 (p) 800.535.6663 (f) 315.685.3361