


Document Description: TSB – RVN 3.4 SOFTWARE RELEASE	Document Number: 20014387 Version: A
	
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.	


<h2>Technical Service Bulletin</h2>			
Product:	RetinaVue Network		Date: 2017-02-01
Subject:	Software Version 3.4 Release		
HW Version(s) Affected:	RV100, RV100-B and RV100-6 with firmware versions 5.00.05 and 5.00.07	SW Version(s) Affected:	RVN 3.0, 3.1, 3.2 or 3.3
Serial Numbers Affected:	N/A	Lot or Date Code Affected:	N/A
Classification:	Informational Only		
Distribution:	<input checked="" type="checkbox"/> Customer Care <input checked="" type="checkbox"/> Product Service <input checked="" type="checkbox"/> Field Service <input type="checkbox"/> ASPs <input type="checkbox"/> Distributors <input type="checkbox"/> Customers <input type="checkbox"/> Company Confidential		
<u>THIS DOCUMENT CONTAINS CONFIDENTIAL INFORMATION</u> It can be distributed or made available to only those individuals, companies, and organizations which have current Welch Allyn Confidentiality Agreements, and have a need to know.			
Training Required:	<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No		
Summary:	Release of RVN 3.4 software		

RetinaVue Network software version 3.4 has been released for production and service. This update applies to customers using version 3.0, 3.1, 3.2 or 3.3. The software changes in version 3.4 include:

1. Application and portal User Accounts secure password policy such as strength and uniqueness requirement, and fixed expiration period.
2. Password resets by sending a temporary password via email.
3. Improving User audit logging traceability.
4. User interface improvements such as contact information input fields to support international formats and exam over-read page displays comments from EMR order if present.

Refer to DIR 60081275 RVN 3.4 Software Release Notes for details of the changes.

Access to the software version 3.4 is through the customer portal.

Document Description: TSB – RVN 3.4 SOFTWARE RELEASE	Document Number: 20014387 Version: A
	
Printed or electronic versions of this document not accessed directly from the designated Welch Allyn Controlled Quality Information System are For Reference Only.	

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
A	N/A	Initial Release	2017-02-01	SC