



The Welch Allyn SureTemp® Plus 690 and SureTemp® Plus 692

- Combines the speed of ear thermometry with the accuracy and reliability of traditional oral thermometry. Takes approximately 4- to 6-second oral, 10-second pediatric axillary, 15-second adult axillary and 10-second rectal temperatures
- Ergonomic design—comfortable to hold, easy to use, and very durable
- Interchangeable, removable probe wells reduce the risk of cross-contamination
- Large backlit LCD display with location identification capability (Model 692 only)
- Includes a 60-second pulse timer for counting pulses and respirations (Model 692 only)
- Your choice of security alarm systems to help prevent theft (Model 692 only)
- Last temperature recall
- Convenient storage for 25 probe covers

WelchAllyn®

Advancing Frontline Care™

THERMOMETERS

US

01692-200	SureTemp Plus Electronic Thermometer (Model 692), Wall Mount, Security System with ID Location Field, and Oral Probe with Oral Probe Well, 3-year limited warranty	01692-MC	SureTemp Plus Electronic Thermometer (Model 692) with Soft Case, Oral Probe and Well, Rectal Probe and Well, 100 PC, 3-year limited warranty
01692-201	SureTemp Plus Electronic Thermometer (Model 692), Wall Mount, Security System with ID Location Field, and Rectal Probe with Rectal Probe Well, 3-year limited warranty	01690-201	SureTemp Plus Electronic Thermometer (Model 690), Rectal Probe with Rectal Probe Well, 2-year limited warranty
01692-300	SureTemp Plus Electronic Thermometer (Model 692), One Per Room, Wall Mount, 9 ft Oral Probe with Oral Probe Well, 3-year limited warranty (must be locked to the wall)	01690-300	SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9 ft Oral Probe with Oral Probe Well, 2-year limited warranty (must be locked to the wall)
01692-301	SureTemp Plus Electronic Thermometer (Model 692), One Per Room, Wall Mount, 9 ft Rectal Probe with Rectal Probe Well, 3-year limited warranty (must be locked to the wall)	01690-301	SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9 ft Rectal Probe with Rectal Probe Well, 2-year limited warranty (must be locked to the wall)
01692-700	SureTemp Plus Electronic Thermometer (Model 692), Wall Mount, 9 ft Oral Probe with Oral Probe Well, 3-year limited warranty, with Rolling Stand and Spare Probe and Well Holder	01690-400	SureTemp Plus Electronic Thermometer (Model 690), Wall Mount, Oral Probe with Oral Probe Well, 2-year limited warranty
01690-200	SureTemp Plus Electronic Thermometer (Model 692), Oral Probe with Oral Probe Well, 2-year limited warranty	01690-700	SureTemp Plus Electronic Thermometer (Model 690), Wall Mount, 9 ft Oral Probe with Oral Probe Well, 2-year limited warranty, with Rolling Stand and Spare Probe and Well Holder

EMEA, ASIA PACIFIC, LATIN AMERICA, CANADA

01692-400	SureTemp Plus Electronic Thermometer (Model 692), Wall Mount, Security System with ID Location Field, and Oral Probe with Oral Probe Well, 3-year limited warranty	01690-410	SureTemp Plus Electronic Thermometer (Model 690), Oral Probe with Oral Probe Well, 2-year limited warranty
01692-401	SureTemp Plus Electronic Thermometer (Model 692), Wall Mount, Security System with ID Location Field, and Rectal Probe with Rectal Probe Well, 3-year limited warranty	01690-401	SureTemp Plus Electronic Thermometer (Model 690), Rectal Probe with Rectal Probe Well, 2-year limited warranty
01692-500	SureTemp Plus Electronic Thermometer (Model 692), One Per Room, Wall Mount, 9 ft Oral Probe with Oral Probe Well, 3-year limited warranty (must be locked to the wall)	01690-500	SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9 ft Oral Probe with Oral Probe Well, 2-year limited warranty (must be locked to the wall)
01692-501	SureTemp Plus Electronic Thermometer (Model 692), One Per Room, Wall Mount, 9 ft Rectal Probe with Rectal Probe Well, 3-year limited warranty (must be locked to the wall)	01690-501	SureTemp Plus Electronic Thermometer (Model 690), One Per Room, Wall Mount, 9 ft Rectal Probe with Rectal Probe Well, 2-year limited warranty (must be locked to the wall)
		01690-501V	SureTemp Plus Veterinary Electronic Thermometer (Model 690), Wall Mount, 9 ft Rectal Probe with Rectal Probe Well, 2-year limited warranty

ACCESSORIES AND PROBE COVERS

WORLDWIDE

05031-101	SureTemp/SureTemp Plus Disposable Probe Covers (case of 1,000 probe covers)	02893-000	Probe & Well Kit 4 ft Oral
05031-105	SureTemp/SureTemp Plus Disposable Probe Covers (case of 5,000 probe covers)	02893-100	Probe & Well Kit 9 ft Oral
05031-125	SureTemp/SureTemp Plus Disposable Probe Covers (case of 1,250 probe covers)	02891-0000	Probe Well, Blue, Model 690/692
05031-150	SureTemp/SureTemp Plus Disposable Probe Covers (case of 1,500 probe covers)	02891-1000	Probe Well, Red, Model 690/692
05031-110	SureTemp/SureTemp Plus Disposable Probe Covers (case of 10,000 probe covers)	08273-000	Thermometer Stand, Model 690/692 with Spare Probe & Well Holder
21326-0000	690/692 Wall Holder	21333-0000	Bracket, Model 600 Stand Upgrade, Model 690/692 with Spare Probe & Well Holder
21330-0000	Holder for Spare Probe & Probe Well, Model 690/692	98024-0000	Extended Warranty—Model 690
01638-000	Calibration Key, Assy, Model 690/692	98025-0000	Extended Warranty—Model 692
02892-000	Probe & Well Kit 4 ft Rectal	01802-110	Welch Allyn/Braun Black Body Calibration Device, 110 V
02892-100	Probe & Well Kit 9 ft Rectal	01802-220	Welch Allyn/Braun Black Body Calibration Device, 220 V

SPECIFICATIONS

Ambient Operating Temperature Range	50.0 °F to 104.0 °F (10.0 °C to 40.0 °C)
Operating Altitude	-557 ft to +16,000 ft (-170 m to +4877 m)
Patient Temperature Range	80.0 °F to 110.0 °F (26.7 °C to 43.3 °C)
Temperature Predict Time*	Oral: Approximately 4-6 secs Adult axillary (18 years and older): Approx. 12-15 secs Pediatric axillary (17 years and younger): Approx. 10-13 secs Rectal: Approx. 10-13 secs
Transport/Storage Temperature	-13 °F to +131 °F (-25 °C to +55 °C)
Calibration Accuracy	+/-0.2 °F (0.1 °C) (Monitor Mode)
Humidity	15% to 95% noncondensing
Dimensions	8.46" x 3.18" x 2.43" (215 mm x 81 mm x 62 mm)
Weight	12.6 ounces (357 grams)
Power	(3) 1.5 VDC AA batteries
Equipment Classification	Class 2 Device, Continuous Operation Not AP or APG equipment
IPXØ	Not protected against the ingress of water

* Predict time depends on probe placement and patient condition. Predict times do not include tissue contact-detection time or time needed to place the probe into the measurement site.

CLASSIFICATIONS

Type of protection against electric shock	Internally powered equipment
Degree of protection against electric shock	Type BF applied part
Degree of protection against ingress of liquids	Not protected
Degree of safety in the presence of flammable anesthetics	Equipment not suitable for use with flammable anesthetics
Mode of operation	Continuous

Welch Allyn, 4341 State Street Road, Skaneateles Falls, NY 13153
(p) 1.888.535.6663 www.welchallyn.com

WelchAllyn[®]

Advancing Frontline Care[™]