

Welch Allyn Ear Care Products



Digital MacroView™ Otoscope

Offers you clear, live and still images of the ear canal and tympanic membrane. Plug in. Power On. Play.

Features and Options

- Improve communication for better patient compliance
- Store digital images in patient records
- Share images with other providers as visual reference
- Excellent training aid for educational purposes
- View small details of tubes and hearing aids
- Special features include: USB cable for connection to your computer; insufflation port for pneumatic otoscopy; adjustable focusing wheel; simple button functions: capture, flip and zoom image; digital imager for megapixel still images; compatible with Welch Allyn 3.5 V power sources

Part Number

23920 Digital MacroView Otoscope with 10 ft USB cable and set of four reusable ear specula



MacroView™ Otoscope

Delivers almost twice the field of view at higher magnification compared to standard otoscopes, making it easier to visualize the auditory canal and tympanic membrane. This instrument is a significant advance in hand-held otoscopy, providing an improved view for enhanced learning and decision making.

Features and Options

- A nearly complete view of the tympanic membrane with approximately twice the field of view and 30% greater magnification than a traditional otoscope
- Better clarity and definition of landmarks
- The ability to adjust focus for variable ear canal length or farsighted eyes
- Tip grip for secure fastening and ejection of ear specula
- Fiber optics produce cool light with no reflections, no obstruction

Part Number

23810 MacroView Otoscope
23820 MacroView Otoscope with Throat Illuminator



Ear Wash System

Welch Allyn's Ear Wash System with Hydrovac Action is a state-of-the-art concept in cerumen removal for patients of all ages.

Features and Options

- Clean—disposable ear tip seals tightly in the ear while outflow is directed into sink, away from you and your patient
- Simple—easy for any healthcare professional to use
- Safe—water temperature and pressure automatically controlled
- Portable—connects to any standard faucet, no batteries or electricity required
- Reimbursable—procedure may qualify for third-party reimbursement

Part Number

29350 Ear Wash System

CPT Codes

69210 Removal of impacted cerumen (separate procedure), one or both ears



MicroTym®3 Portable Tympanometric Instrument

Portable, reliable, cost-effective tympanometry for patients of all ages.

Features and Options

- 1000Hz quick interpretation mode for testing of infants < 6 months
- Enhanced, field upgradeable software
- Use for the diagnosis of otitis media and other middle ear conditions such as otosclerosis, cholesteatoma and patent tympanostomy tube
- Test in one second, print out results in five seconds
- LCD preview screen for reliable results

Part Number

93650 MicroTym3 Set with Charging Stand and Printer
23650 MicroTym3 Set Handle Only (includes handle, battery and tip box with 2 tips of each size)

CPT Codes

92567 Tympanometry (impedance testing)

WelchAllyn®

Advancing Frontline Care™



TM 286™ AutoTymp® Tympanometer

More data for a more complete diagnosis.

Features and Options

- Three-in-one tympanometry, acoustic reflex and audiometry testing
- Provides multiple frequencies (125Hz to 8000Hz)
- Wide intensity range
- USB for connectivity
- Portable design

Part Number

28600 TM 286 AutoTymp

CPT Codes

- 92551 Screening test, pure tone, air only
- 92552 Pure tone audiometry (threshold); air only
- 92567 Tympanometry (impedance testing)
- 92568 Acoustic reflex testing; threshold



AM 282™ Manual Audiometer

Threshold audiometry to quantify hearing loss.

Features and Options

- Quantifies specific dB and frequency hearing loss with precision testing
- Provides multiple frequencies (125Hz to 8000Hz)
- Wide intensity range
- Compact, light-weight

Part Number

28200 AM 282 Manual Audiometer (with case)

CPT Codes

- 92551 Screening test, pure tone, air only
- 92552 Pure tone audiometry (threshold); air only



AudioScope® 3 Screening Audiometer

View the tympanic membrane and ear canal while screening for hearing loss, quickly and accurately.

Features and Options

- Halogen light and fiber optics enable brilliant view of tympanic membrane and ear canal before and during screening
- Choice of three screening levels—20 dB, 25 dB HL, and 40 dB HL—at 500, 1000, 2000 and 4000 Hz for all ages
- Built-in 1000 Hz pretone lets patients “practice” listening first

Part Number

92680 AudioScope 3 Set (AudioScope 3, set of AudioSpec® Specula, 110-130 V Charging Transformer, Charging Stand and Recording Forms)

CPT Codes

- 92551 Screening test, pure tone, air only



OAE Hearing Screener

Precisely and objectively screen for hearing loss in newborns, infants, toddlers, preschool and school-age children in just seconds.

Features and Options

- Requires minimal patient cooperation
- Fast, accurate results (about 10 seconds per ear)
- Objective testing—no patient response required
- Clear results that don't require caregiver interpretation
- Portable design
- 10-test memory storage in device
- Simple printer interface for easy patient data documentation
- Download patient information into PC ASCII file
- Optional power supply: rechargeable lithium battery or AC power
- One-year warranty, extended warranty available
- Fits a variety of ear tip sizes

Part Number

29400 Welch Allyn OAE Hearing Screener Set, including Printer

CPT Codes

- 92587 Evoked otoacoustic emissions; limited (single stimulus level), either transient or distortion products

*CPT is a registered trademark of the American Medical Association

CPT Codes © 2008 American Medical Association. All rights reserved. CPT is a trademark of the AMA. No fee schedules, basic units, relative values or related listings are included in CPT. The AMA assumes no liability for the data contained herein.

Applicable FARS/DFARS restrictions apply to Government use.

For reference only. Information does not constitute a guarantee of coverage or payment.

WelchAllyn®

Advancing Frontline Care™