

23920 Digital MacroView Magnification

Magnification of the Digital MacroView is dependent on the size of the viewing monitor and the program that is being used to view the image (MS Windows Explorer, Welch Allyn Video Viewer, etc).

The values presented below assume the use of a 14.1 inch laptop display, Windows Explorer and the Digital MacroView set at the nominal focus position.

- Nominal Magnification = 16.6X
- Nominal magnification with zoom = 33.2X

Under the same conditions:

- Focused short, mag = 24X or 48X zoomed
- Focused long, mag = 11X or 22X zoomed

Note: in the zoomed state the image is approximately ¼ of the standard (non-zoomed) view