

Computer Specifications for RetinaVue™ Network Client:

The following are the minimum requirements for any workstation that will run the Welch Allyn RetinaVue Network Client Application.

Supported Operating Systems and Memory*:

- Windows XP Pro Service Pack 3 (32-bit) with at least 1 GB of RAM
- Windows 7 Pro or Ultimate (64-bit) with at least 4 GB of RAM
- Windows 8 Pro (64-bit) with at least 4 GB of RAM

* More memory may be required depending on requirements of other applications installed on the computer.

Processor: 1.5 GHz or higher

Disk Space: 150 MB or more

Screen Resolution: 1280 x 720 or higher

Universal Serial Bus: Must have at least 1 available USB port

Ethernet Adapter: Must have an Ethernet port to accept a network cable (RJ-45 connector)

Web browser:** Microsoft Internet Explorer 8 or higher

** To ensure compatibility as browser technologies evolve, Welch Allyn may change the supported browsers set forth above from time to time as older versions become obsolete and newer versions become available.