



Welch Allyn Inc.

Medical Division

4341 State Street Road, Box 220  
Skaneateles Falls, New York 13153-0220 USA  
Tel: (315) 685-4445 FAX: (315) 685-4653

## Technical Specification Bulletin

*This Bulletin references, and adds clarification to Product specification for:*

### **Otosopes ... optical specifications**

Model	Description	Lens Part Nbr.	Diopters	Focal Length	Mag.	Window
20000 / 25020	3.5V Otoscope	200057	8.8	114 mm 4.5 in.	2.2 x	Black
20100/ 20200	Standard Otoscope	201016	10	100 mm 4 in.	2.5 x	Black
20160 / 20260	Vet. Otoscope	235025	8	125 mm 5 in.	2.0 x	Green
20151 / 20251	12 Diopter Otoscope	201016-1	12	83 mm 3.5 in.	3.0 x	Grey
21110	2.5V Pocketscope	211005-1	9.8	102 mm 4 in.	2.5 x	Clear
21600 / 21700	Operating Otoscope	209009	7.9	127 mm 5 in.	2.0 x	Black
21660 / 21760	Vet Operating Otoscope	209021	6.5	152 mm 6 in.	1.7 x	Green
20202	Consulting Otoscope	202114	6.5	153 mm 6 in.	1.7 x	Black
20262	Vet Consulting Otoscope	202116	5.6	177 mm 7 in.	1.4 x	Green
23540	Mag Lens Otoscope	328075-2	2	500 mm 19.7 in.	0.5 x	Black
24031	2.5 v Otoscope	200057	8.8	114 mm 4.5 in.	2.2 x	Black
29090	Otoscope	290005	9.8	102 mm 4 in.	2.5 x	Clear

Construction Materials:

ABS and Acrylic Plastics, Stainless Steel, Chrome Plated Brass, Glass.