

# Welch Allyn


## *Tycos*

### *Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog*

PN 95P507 Rev B.

© 2001 Welch Allyn Inc.

All Rights Reserved



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Revision History

Revision	Date	Author	Description
A	1/24/01	JDB	Initial Release of Catalog
B	10/12/01	JDB	Removed the columns that listed the prices of the parts.

# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## *Table Of Contents*

Section 1: Classic Pocket Aneroid Replacement Parts .....	4
Section 2: Classic Hand Aneroid Replacement Parts .....	5
Section 3: 767 Wall/Mobile Aneroid Replacement Parts .....	6
Section 4: Classic Wall/Mobile Aneroid Replacement Parts.....	7
Section 5: Economy Aneroid Replacement Parts .....	8
Section 6: Silver Ring Hand Aneroid Replacement Parts.....	9
Section 7: Pro Check TR2 Hand Aneroid Replacement Parts .....	10
Section 8: The Elite Stethoscope Replacement Parts.....	11
Section 9: Harvey DLX Stethoscopes Replacement Parts .....	13
Section 10: Original Harvey Stethoscopes Replacement Parts .....	14
Section 11: Spectrum Stethoscopes Replacement Parts .....	15
Section 12: Medallion Stethoscopes Replacement Parts .....	17
Section 13: Meditron Electronic Stethoscope Replacement Parts .....	18

# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 1: Classic Pocket Aneroid Replacement Parts

The following are the list of replacement parts for the Pocket Aneroid Gauge that can be purchased through Welch Allyn Inc.

<u>Part Number</u>	<u>Description</u>
51P52	Clip Assembly
2P1623	Case
32P180	Lock nut
34P760A	Crystal
3P831	Bezel
35P804	Screw cap (NS)
61S681	Diaphragm Assembly
63S401	Pointer
93S617	Movement
9P1799	Screw, Movement
11P363	Spring
7P689	Pin for Clip



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 2: Classic Hand Aneroid Replacement Parts

The following are the list of replacement parts for the Classic Hand Aneroid Gauge that can be purchased through Welch Allyn Inc.

<u>Part Number</u>	<u>Description</u>
12P368	Lure Lock Fitting
9P2651	Screw, Cover
9S132	Thumb Screw
42P252	Valve Body Assembly
5090-30	Screen Assy/Bulb Adapter
12P490	Bulb Adapter
2P1627A	Case
5086-01	Bulb
2P1628	Cover, Hand Gauge
3P831	Bezel
34P760A	Crystal
63S401	Pointer
9P1799	Screw, Movement
93S617	Movement
61S694	Diaphragm Assembly
32P527	Lock Nut



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 3: 767 Wall/Mobile Aneroid Replacement Parts

The following are the list of replacement parts for the 767 Wall/Mobile Aneroid Gauge that can be purchased through Welch Allyn Inc.

### 767 Model

<u>Part Number</u>	<u>Description</u>
3P844	Bezel
34P768	Crystal
63S412	Pointer
93S617A	Movement
61S681	Diaphragm Assembly
2P1138	Mount
32P727	Nut
2P1631	Case
32P763	Lock nut, Case - <i>Not pictured</i>
8P1064	Bushing - <i>Not pictured</i>
96P135	Basket - <i>Not pictured</i>
9P2841	Dial Mounting Screw



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 4: Classic Wall/Mobile Aneroid Replacement Parts

The following are the list of replacement parts for the Classic Wall/Mobile Aneroid Gauge that can be purchased through Welch Allyn Inc.

### Classic Model

Part Number	Description
2P1138	Mount
9P1804	Spacer
61S681	Diaphragm Assembly
32P727	Jam Nut
14P1718	Bracket
9P2455	Bracket Screw
8P949	Bracket Bushing
32P763	Bracket Nut
93S617A	Movement
34P417	Crystal
9P2841	Dial Mounting Screw
63S298	Pointer, USA
63S402	Pointer, Japan - <i>Not pictured</i>
2P1503	Case (Gray)
2P1124	Case (Beige) - <i>Not pictured</i>
9P2465	Screw (Case)
22P830	Washer (Case)



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 5: Economy Aneroid Replacement Parts

The following are the list of replacement parts for the Economy Aneroid Gauge that can be purchased through Welch Allyn Inc.

<u>Part Number</u>	<u>Description</u>
2P1437	Cover for Economy Aneroid
2P1626C	Case for Economy Aneroid
34P640	Crystal
93S618	Movement
2P1130D	Clip
32P727	Jam Nut
63S401	Pointer
9P2563	Screw (Case)
61S693	Capsule
9P1799	Screw, Movement





# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 6: Silver Ring Hand Aneroid Replacement Parts

The following are the list of replacement parts for the Silver Ring Hand Aneroid Gauge that can be purchased through Welch Allyn Inc.

Part Number	Description
12P368	Adapter
75P240	O-Ring
2P1632	Housing
33P131	Movement Cup
70P450	Button
7P716	Pivot Pin
93S720	Movement
22P1455	Washer Silicon
34P434	Crystal
39P505	Trim Ring
52P1848	Dial
18P189	Bulb Collar
23P428	Bulb Adapter
97S581	Bulb
11P1240	Spring
42S242	Valve Body
8S61	Valve Assembly
5090-30	Screen Assy/Bulb Adapter
3P847	Bezel
63S410	Pointer



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 7: Pro Check TR2 Hand Aneroid Replacement Parts

The following are the list of replacement parts for the Pro Check TR2 Hand Aneroid Gauge that can be purchased through Welch Allyn Inc.

<u>Part Number</u>	<u>Description</u>
34P761	Crystal
63S416	Pointer
52P1845	Dial
93S720	Movement
22P1460	Washer
97S583	Bulb
18P190	Bulb Ring
70P446	Trigger
8S61	Valve Assy
11P1240	Spring - <i>Not pictured</i>
52-09-950	<i>Housing (order from S+K)</i>
15-40-000	<i>Valve Body (order from S+K)</i>



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 8: The Elite Stethoscope Replacement Parts

The following are a list of spare parts for the Elite Stethoscope

Part Number	Description
85S270	Tubing, Black, 22"
85S269	Tubing, Black, 25"
85S269B	Tubing, Burgundy, 25"
85S269N	Tubing, Navy, 25"
85S269FG	Tubing, Forrest, 25"
51S174	Binaural/Spring Assy, Black
51S174B	Binaural/Spring Assy, Burgundy
51S174N	Binaural/Spring Assy, Navy
51S174FG	Binaural/Spring Assy, Forrest Green
42P389	Chest piece Pediatric,
42P390R	Chest piece, Adult,
3P837	Non-chill Rim, Bell, Pediatric, Black
3P837G	Non-chill Rim, Bell, Pediatric, Gray
3P836	Non-chill Rim, Bell, Adult, Black
3P836G	Non-chill Rim, Bell, Adult, Gray
3P850	Non-chill Rim, Diaphragm, Pediatric, Black
3P850G	Non-chill Rim, Diaphragm, Pediatric, Gray
3P849	Non-chill Rim, Diaphragm, Adult, Black
3P849G	Non-chill Rim, Diaphragm, Adult, Gray
61P373	Diaphragm, Pediatric, Gray
61P374	Diaphragm, Pediatric, Black
61P375	Diaphragm, Adult
61P383K	Diaphragm, WA, Black Holographic
61P383N	Diaphragm, WA, Navy Holographic
61P383B	Diaphragm, WA, Burgundy Holographic
61P383FG	Diaphragm, WA, F. Green Holographic

Part Number	Description
45S100	Name Tag
45p3239	Plate for Name Tag
5079-113	Hard Eartips ( x-large, pair)
5079-66	Hard Eartips ( large, pair)
5079-233	Comfort Sealing Eartips, Gray
5079-336	Comfort Sealing Eartips, Black
5079-337	Comfort Sealing Eartips, Navy
5079-	Comfort Sealing Eartips, Burgundy
5079-	Comfort Sealing Eartips, Forest Green

# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 9: Harvey DLX Stethoscopes Replacement Parts

The following are a list of spare parts for the Harvey DLX Stethoscopes

<u>Part Number</u>	<u>Description</u>
45S100	Name Tag
45P3239	Plate for Name Tag
85S269	Tubing 25", Black
85S269N	Tubing 25", Navy
85S269B	Tubing 25", Burgundy
85S269FG	Tubing 25", Forrest Green
51S174	Binaural/Spring Assembly, Black
51S174N	Binaural/Spring Assembly, Navy
51S174B	Binaural/Spring Assembly, Burgundy
51S174FG	Binaural/Spring Assembly, F. Green
35P1339	Hard Eartips, X-large, Pair
35P779	Hard Eartips, Large, Pair
5079-233	Comfort Sealing Eartips, Gray
5079-336	Comfort Sealing Eartips, Black
5079-337	Comfort Sealing Eartips, Navy
5079-338	Comfort Sealing Eartips, Burgundy
5079-339	Comfort Sealing Eartips, Forrest Green
61P386K	Diaphragm, Harvey Holographic, Blk
61P386N	Diaphragm, Harvey Holographic, Navy
61P386B	Diaphragm, Harvey Holographic, Burg
61P386FG	Diaphragm, Harvey Holographic, F. Green



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 10: Original Harvey Stethoscopes Replacement Parts

The following are a list of spare parts for the Harvey Original stethoscopes.

Part Number	Description
5079-20	Harvey Accessory Kit
47P1108	Tubing, (two tubes joined for use with two-piece spindle)
47P-1117	Tubing, (two tubes joined for use with one-piece spindle)
47P1124	Tubing, (two pieces, for use with two-piece spindle) *requires two each
61S541	Flat Adult Chest Piece
61S540	Corrugated Adult Chest Piece
42S195	Flat Pediatric Chest Piece
61P363	Diaphragm, Pediatric
42P175	Bell Chest Piece
61P301	Flat Diaphragm only
61P362	Corrugated Diaphragm only
61P363	Pediatric Diaphragm only
51S16	Binaurals
45P2967	Name Clip
35P1214	Eartips (mushroom)
5079-233	Comfort Sealing Eartips
5079-111	Soft Eartips (medium, pair)
5079-110	Soft Eartips (large, pair)
35P1216	Hard Eartips (small, pair)
35P779	Hard Eartips (large, pair)
35P1339	Hard Eartips (x-large, pair)



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 11: Professional Stethoscopes Replacement Parts (continued next pg.)

The following are a list of spare parts for the Spectrum Single Head Stethoscopes

### Non-Chill Rings

Part number	Description
3P841K	Adult Bell, Black
3P841G	Adult Bell, Gray
3P842K	Pediatric Bell, Black
3P842G	Pediatric Bell, Gray
3P849	Adult Diaphragm, Black
3P849G	Adult Diaphragm, Gray
3P850	Pediatric Diaphragm, Black
3P850G	Pediatric Diaphragm, Gray

### Eartips (requires two each tip and insert per stethoscope)

Tip	Description
35P1346	Medium Plastic, pair
35P1347	Large Plastic, pair
5079-170	Comfort Sealing, Gray
5079-366	Comfort Sealing, Black
5079-367	Comfort Sealing, Navy
5079-368	Comfort Sealing, Burgundy
5079-369	Comfort Sealing, Forest Green

### Diaphragm Discs

Part number	Description
61P375	Adult
61P374	Pediatric
61P383K	Diaphragm, Holographic, Adult, Blk
61P383N	Diaphragm, Holographic, Adult, Navy
61P383B	Diaphragm, Holographic, Adult, Burgundy
61P383FG	Diaphragm, Holographic, Adult, F. Green



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 11: Professional Stethoscopes (cont.)

The following is a continuation of the spare parts list for the Spectrum Single/Double Head Stethoscopes

### Double Head Chestpiece Assemblies

Part Number	Description
45S251R	Adult, Chest Piece
45S252	Pediatric, Chest Piece
3P849G	Ring, Adult, Gray
3P849	Ring, Adult, Black
3P850G	Ring, Pediatric, Gray
3P850	Ring, Pediatric, Black
3P841G	Bell Ring, Adult Gray
3P841K	Bell Ring, Adult Black
3P842G	Bell Ring, Pediatric Gray
3P842K	Bell Ring, Pediatric Black
61P375	Diaphragm Adult (single and double)
61P374	Diaphragm, Ped (single and double)

### Single Head Chestpiece Assemblies

42S245R	Adult, Chest Piece
42S246	Pediatric Chest Piece
3P849	Adult, Black Ring
3P849G	Adult, Gray Ring
3P850	Pediatric, Black Ring
3P850G	Pediatric, Gray Ring

### Y-Tube/Binaural Assemblies

85P233K	Black 23" (58.4 cm)
85P233G	Gray 23" (58.4 cm)
85P233N	Navy 23" 58.4 cm)
85P233B	Burgundy 23" (58.4 cm)
85P233FG	Forrest Green 23" (58.4 cm)
85P233BL	Blue 23" (58.4 cm)
85P233K	Black 28" (71.1 cm)
85P233G	Gray 28" (71.1 cm)
85P233N	Navy 28" (71.1 cm)
85P233B	Burgundy 28" (71.1 cm)
85P233FG	Forrest Green 28" (71.1 cm)
85P233BL	Blue 28" (71.1 cm)
51S176	Binaural, Spectrum

Note: 23" Tubing is made form the 28" tubing. Thus the same part number.



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 12: Medallion Stethoscopes Replacement Parts

The following are a list of spare parts for the Medallion Gift Stethoscopes

<u>Part Number</u>	<u>Description</u>
42S194G	Adult Chestpiece Assembly
42S195G	Pediatric Chestpiece Assembly
42P175G	Bell Chestpiece Assembly
5079-353	Diaphragm, Adult, Holographic, 4 colors



# Welch Allyn Aneroid Sphygmomanometer and Stethoscope Spare Parts Catalog

## Section 13: Meditron Electronic Stethoscope Replacement Parts

The following are a list of spare parts for the Meditron Electronic Stethoscope.

<u>Part Number</u>	<u>Description</u>
39P507	Ring, Speaker housing
167S105	Speaker
35P1390	Front Speaker housing
94P2667	Front Foam, Speaker Housing
94P2668	Rear Foam, Speaker housing
114P85	Battery (2 ea.)
3P853	Battery Cover
35P1392	Ear Tips
66P821	Cables for Distributor
5079-414	Headset

