

Stethoscope ID Matrix				
Electronic Stethoscopes	Identifiers		Acoustic Stethoscopes	Identifiers
Master Elite	Obsolete 12/31/2011			
5079-405	Picture of Heart and Lung		Harvey DLX	
5079-405V	Picture of Heart and Lung Requires 1 cylindrical battery 5079-425		5079-321, 5079-322, 5079-323, 5079-324, 5079-325, 5079-326 5079-327, 5079-328, 5079-325P	No non chill rims, single tube from patient end of scope, coming up to a "V" then a hard plastic "U"
			Harvey Elite	
Master Elite Plus			5079-122, 5079-125, 5079-270, 5079-271, 5079-284, 5079-125P	Non chill rims black, single tube from the patient end of the scope, coming up to a "V" then a hard plastic "U". May have a date code (Z97) indicating the year of manufacture. This would be found on the end of the tubing near the bell and diaphragm.
5079-400	Obsolete 1/31/2008			
5079-400V	Obsolete 1/31/2008		Professional Adult aka Spectrum	
	Required 2 flat disc batteries 5079-410		5079-135, 5079-137, 5079-289, 5079-139, 5079-285	Non chill rims black or gray, single tube coming up to a "U" in the binaurals. Diaphragm 1+7/8 across, bell 1+3/8
TrimLine			Professional Pediatric aka Spectrum	
Echo Electronic	Not Supported by Welch Allyn		5079-145, 5079-147, 5079-149, 5079-287	Non chill rims black or gray, single tube coming up to a "U" in the binaurals.
Symphony Tunable	Not Supported by Welch Allyn			
Puretone Cardiology	Not Supported by Welch Allyn		Lightweight Double Head	
Puretone Traditional	Not Supported by Welch Allyn		5079-73, 5079-76, 5079-74, 5079-75	Single tube from the patient end of the scope, coming up to a "U" with a chrome plated spring in the binaurals. No parts or accessories.
Symphony Traditional	Not Supported by Welch Allyn			
Economy	Not Supported by Welch Allyn		Disposable Uniscope	
			17462P blue, 17461 yellow	All plastic. No parts or accessories, MRI compatible. P = Pediatric
			Obsolete Stethoscopes	Identifiers
			Original Harvey	
			5079-28,5079-27	Two tubes separate or fused together from the patient end of the scope splitting near the binaural with a plastic coated spring. No non chill rims.
			Medallion Gift (2003)	Blue box
			Anesthesiologist Pediatric (2003)	Single head looks the same as the Professional (except it is a single head not a double head) Tubing separates in the middle.
			Anesthesiologist Adult (1/31/2008)	Single head looks the same as the Professional (except it is a single head not a double head) Tubing separates in the middle.