


FlexiPort Lot Code-2013 on

Pre 2013, the Lot Code was on back of cuff, above the velcro - Now it is on front of cuff, above the people.

- The FlexiPort lot code is imprinted in the top right corner of the cuff



- The lot code is in the format YYYY-MM-DD-XX
(ex. Lot 2013-10-15-12 was made on October 15, 2013, the last digits (12) identify the machine that they were made on).



FlexiPort Lot Code

- The FlexiPort lot code is imprinted in the top right corner of the cuff



- The lot code is in the format YYYY-MM-DD-XX
(ex. Lot 2013-08-01-09 was made on August 01, 2013, the last digits are the machine that they were made on).

