

Hear the Difference

from the leader in diagnostic care.

FOR VETERINARIANS

Welch Allyn Veterinary Stethoscopes

For over a decade, we have been focused on providing veterinarians with exceptional performance and value from every Welch Allyn device. Our stethoscopes are no exception.



Harvey™ DLX™ Stethoscopes

Trumpet brass construction translates to hearing more of the important cardiac, pulmonary, vascular and abdominal sounds without relying on difficult “variable pressure” techniques.

- Available with double- or triple-head rotatable design
- Chrome-plated brass chestpiece and binaurals
- Rotatable binaurals and interchangeable comfort-sealing ear tips
- Includes pediatric diaphragm for use with small animals
- 28" tubing length



5079-315 Accessory Kit

(Includes one pediatric chestpiece, one flat diaphragm, one corrugated diaphragm, one pair soft eartips, one pair hard eartips)



Harvey™ Elite® Stethoscope

Superior comfort and a stainless steel head that transmits sounds crisply and clearly.

- Double head (bell and diaphragm)
- Rotatable binaurals and interchangeable comfort-sealing ear tips
- Pediatric diaphragm available
- 28" tubing length



5079-232 Accessory Kit

(Includes one pediatric chestpiece, one flat diaphragm, one corrugated diaphragm, one pair soft eartips, one pair hard eartips)

Product Description	HARVEY DLX	HARVEY ELITE
Product Type	Double/Triple Head	Double Head
Head Diameter	1.9" (48.3 mm)	1.9" (48.3 mm)
Overall Length	28" (711 mm)	28" (711 mm)
Tubing Length	16.5" (419 mm)	16.5" (419 mm)
Latex Free	Yes	Yes
Chestpiece Material	Brass	Stainless Steel
Diaphragm Material	Polycarbonate	Polycarbonate
Diaphragm Style	Flat or Corrugated	Flat or Corrugated
Type of Tubing	Dual Lumen	Dual Lumen
Anti-Chill Ring	No	Yes
Warranty	10-year	10-year
Product Numbers	<ul style="list-style-type: none"> ● 5079-321 Triple, Black ● 5079-323 Triple, Navy ● 5079-325 Double, Black ● 5079-327 Double, Navy ● 5079-325P Pediatric, Double, Black ● 5079-322 Triple, Burgundy ● 5079-324 Triple, Forest Green ● 5079-326 Double, Burgundy ● 5079-328 Double, Forest Green 	<ul style="list-style-type: none"> ● 5079-125 Elite, Black ● 5079-270 Elite, Burgundy ● 5079-271 Elite, Navy ● 5079-284 Elite, Forest Green ● 5079-125P Pediatric, Elite, Black

To learn more about Welch Allyn Stethoscopes contact your local distributor rep, or call 1.800.535.6663.

4341 State Street Road, Skaneateles Falls, NY 13153 USA (p) 800.535.6663 (f) 315.685.3361

© 2012 Welch Allyn MC9131

www.welchallyn.com/dvm

WelchAllyn®

Advancing Frontline Care™