



Welch Allyn
Connex® Spot Monitor

Welch Allyn Connex[®] Spot Monitor

Accurate, automatic vitals sent to your Electronic Medical Records (EMR) from the point of care for improved patient safety and clinical decisions.

The monitor to turn to for routine vital signs acquisition, interval monitoring, and custom scoring. Plus EMR integration to give clinicians immediate access to the information they need—where and when they need it.

Improve patient outcomes and satisfaction

- > Reduce transcription errors and omissions by automatically transferring patient vital signs and other data directly into your EMR
- > Customisable fields to fully document at the point of care, eliminating the step of also entering some data in computers
- > Compatible with a wide range of cleaning products to reduce the risk of cross-contamination
- > Industry-leading parameter technologies to ensure accuracy and performance across a diverse range of patient populations

Accessory Power Management Stand with WhisperDrive™ technology

- > Exclusive power management features to ensure performance even if the prior shift forgets to recharge the device
- > WhisperDrive technology provides smooth, quiet maneuverability to help minimise patient disturbance
- > LED lighting to help clinicians work in dark rooms without disturbing patients
- > Small work surface to improve ergonomics and ensure there's always a convenient place to document notes
- > Organisation of accessories to help prevent tangling and sensor damage
- > Cleaning supply storage to encourage proper device cleaning

Upgradeable design to help protect your investment

Designed for the ward environment with an easy-to-use, configurable touchscreen display. And future-proof your purchase with these upgrade options designed to help keep patients safe and improve outcomes:

- > Masimo[®] SpO₂
- > Nellcor[®] OxiMax™ SpO₂
- > Nonin[®] SpO₂
- > Braun PRO 6000 thermometry
- > Barcode scanning
- > Compatible with a wide range of wireless networks and security protocols
- > Access to remote software upgrades and USB accessories





Improve staff satisfaction and productivity

- > Designed by clinicians to increase efficiency and staff satisfaction
- > Intuitive touchscreen user interface to simplify training and operation
- > On-screen guides to prompt inexperienced users

Secure patient information

- > Optional clinician sign-on screen restricts access to protected health information
- > Optional features, like abbreviations, help protect patient-identifiable information during clinical use
- > Automated deletion of patient data after sending to EMR



Early Warning Scores

- > Captures custom patient scores at the bedside to allow clinicians to document and track signs of patient deterioration
- > Use upto three different scoring systems, in line with your hospital's guidelines—such as MEWS, PEWS, etc.
- > Displays facility's own recommendations based on patient scoring results to help guide clinicians through next steps and to make informed care decisions quickly



Partners in CareSM Services for Connex Spot Monitor

Welch Allyn offers a full suite of services to help protect your organisation's investment and ensure clinicians have what they need to make informed patient care decisions.

- > **Unique Enterprise Deployment Solution**—to ensure your monitors arrive configured to your specification right out of the box
- > **Biomedical Service and Repair Training**—to empower staff to maintain and repair devices in-house
- > **Support Programs**—cover repair parts and replacement accessories for efficient, cost-effective device support
- > **Service Dashboard**—enables proactive monitoring of connected devices, configuration creation, and remote deployment

Configurations Available

CONNEX SPOT MONITOR	7100	7400	7500
Parameters			
SureBP® NIBP	●	●	●
SureTemp® Plus Thermometry	○	○ ↑	○ ↑
Braun PRO 6000 Ear Thermometry	○	○ ↑	○ ↑
Masimo or Nellcor OxiMax® SpO ₂		●	●
Nonin SpO ₂	○		○
Communications			
WiFi		↑	●
USB, Ethernet	●	●	●
Accessories			
Classic Stand	○	○	○
APM Stand	○	○	○
GCX Wall Channel	○	○	○
Barcode Scanner	○	○	○

● Included
 ○ Optional
 ↑ Upgradeable

ACCESSORIES

7000-APM	Connex Spot Accessory Power Management Stand
7000-MS3	Connex Spot Classic Mobile Stand
7000-GCX	Connex Spot GCX VESA Wall Channel
6000-NC	VSM 6000 Nurse Call Cable
7000-916HS	HS1 2D Barcode Scanner
106275	Connectivity Accessory Kit USB Cable
6000-50	VSM 6000 USB Configuration Memory Stick
7000-PS	Power Supply (required if not ordering a stand or wall mount)
BATT22	CSM Lithium-Ion Battery
BATT99	APM Stand Lithium-Ion Battery
4500-35	Blood Pressure Hose with FlexiPort (10 ft)
7000-33	Neonatal Blood Pressure Hose (10 ft)
LNCS-DCIP	Masimo Pediatric Reusable Finger Sensor
D-YS	Nellcor Dura-Y® Sensor
D-YSPD	Nellcor PediCheck® Pediatric Finger Sensor (Requires D-YS)
2360-010	Nonin Reusable Pediatric Sensor (2 m)
S1-CSM-5	CSM, Comprehensive Partner Program, 5 years
S2-CSM-5	CSM, Biomed Partner Program, 5 years
S4-CSM-5	CSM, Warranty Extension, 5 years

One- and two-year service plans also available.

Welch Allyn Middle East
 Dubai Healthcare City, Unit #407, 4th Floor Building 49
 Dubai, United Arab Emirates - P.O. Box 89106
 (p) +971 4 3623496 (f) +971 4 43686547
 Email: info.middleeast@welchallyn.com

WWW.WELCHALLYN.COM

WelchAllyn®
 Advancing Frontline Care™