



Welch Allyn CP 50™ Electrocardiograph Part Number Configurations



CP50 - _ - _ -

Model Type	Patient Cable Labeling & Electrode Type	Language	Power Cords
A - Analysis	1 - AHA, Disposable	EN - English	1 - US
P - Plus	2 - IEC, Disposable	FR - French	2 - European
Combinations • CP50-XXX • CP50 A -XXX • CP50 P -XXX • CP50 AP -XXX	3 - AHA, Reusable	DE - German	3 - Israel
	4 - IEC, Reusable	ES - Spanish	4 - UK
		IT - Italian	5 - Switzerland
		NO - Norwegian	6 - Australia
		NL - Dutch	7 - China
		SV - Swedish	8 - Denmark
		DA - Danish	9 - South Africa
		PT - Portuguese	
		ZH - Chinese	
		PL - Polish	
	RU - Russian		

Example: **CP50AP-2ZH7**

Language abbreviations match ISO standard.

Accessories			
	Description (Part Number)		
Included (Reference part numbers to reorder)	100 Disposable Tab Electrodes (45008-0000), and 10 Adaptors (715006)	OR	6 Suction Cups (714730), and 4 Limb Clamps (714731)
	AHA Patient Cable (714507)	OR	IEC Patient Cable (714508)
	Z-fold Paper (406021)		
	Dust Cover (713942)		
	Battery (406812)		
Optional	Roll Stand (406815)		
	Carrying Case (715107)		
	CP 50 Connectivity Kit (406814)		